

# ISO 16000-6:2004-03 (E)

## Indoor air - Part 6: Determination of volatile organic compounds in indoor and test chamber air by active sampling on Tenax TA sorbent, thermal desorption and gas chromatography using MS/FID

---

<b>Contents</b>		<b>Page</b>
Foreword .....		iv
Introduction .....		v
1	Scope .....	1
2	Normative references .....	1
3	Terms and definitions .....	2
4	Principle .....	2
5	Reagents and materials .....	3
6	Apparatus .....	4
7	Conditioning and storage of sorbent tubes .....	6
7.1	Conditioning .....	6
7.2	Storage of conditioned sorbent tubes before sampling .....	6
8	Calibration of pump .....	6
9	Sampling .....	6
9.1	Indoor air sampling .....	6
9.2	Test chamber air sampling .....	7
9.3	Sampling volumes .....	7
9.4	Storage of loaded samples .....	7
9.5	Field blanks .....	7
10	Analysis .....	8
10.1	General .....	8
10.2	Thermal desorption .....	8
10.3	Temperature programme .....	8
10.4	Analysis of the samples .....	9
11	Identification of single VOCs .....	9
12	Concentration of analytes in the sampled air .....	9
12.1	General .....	9
12.2	Volatile organic compounds (VOCs) .....	9
12.3	Total volatile organic compounds (TVOCs) .....	10
12.4	VVOC and SVOC compounds observed outside the TVOC range .....	11
13	Performance characteristics .....	11
14	Test report .....	12
15	Quality control .....	13

<b>Annex A (informative) Examples of compounds detected in indoor air and from building products in test chambers .....</b>	<b>14</b>
<b>Annex B (informative) Safe sampling volumes (SSV) for selected organic vapours sampled on Tenax TA .....</b>	<b>20</b>
<b>Annex C (informative) Storage recovery of solvents on Tenax TA® sorbent tubes .....</b>	<b>22</b>
<b>Bibliography .....</b>	<b>24</b>