

# ISO 13099-2:2025-08 (E)

## Colloidal systems - Methods for zeta-potential determination - Part 2: Optical methods

---

<b>Contents</b>		<b>Page</b>
Foreword .....		iv
Introduction .....		v
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Normative references .....</b>	<b>1</b>
<b>3</b>	<b>Terms, definitions and symbols .....</b>	<b>1</b>
<b>3.1</b>	<b>Terms and definitions .....</b>	<b>1</b>
<b>3.2</b>	<b>Symbols .....</b>	<b>3</b>
<b>4</b>	<b>Principles .....</b>	<b>3</b>
<b>5</b>	<b>Microscopic methods .....</b>	<b>4</b>
<b>6</b>	<b>Electrophoretic light-scattering (ELS) method .....</b>	<b>5</b>
<b>6.1</b>	<b>General .....</b>	<b>5</b>
<b>6.2</b>	<b>Cell design .....</b>	<b>5</b>
<b>6.3</b>	<b>Reference beam optical arrangement .....</b>	<b>6</b>
<b>6.4</b>	<b>Signal processing .....</b>	<b>8</b>
<b>6.4.1</b>	<b>Spectrum analysis .....</b>	<b>8</b>
<b>6.4.2</b>	<b>Autocorrelation function .....</b>	<b>8</b>
<b>6.4.3</b>	<b>Phase analysis light scattering (PALS) .....</b>	<b>9</b>
<b>6.4.4</b>	<b>Modulated Brownian motion power spectrum method .....</b>	<b>9</b>
<b>6.5</b>	<b>Determination of electrophoretic mobility .....</b>	<b>10</b>
<b>7</b>	<b>Calculation of zeta-potential .....</b>	<b>10</b>
<b>8</b>	<b>Operational procedures .....</b>	<b>11</b>
<b>8.1</b>	<b>Requirements .....</b>	<b>11</b>
<b>8.1.1</b>	<b>Instrument location .....</b>	<b>11</b>
<b>8.1.2</b>	<b>Dispersion liquids .....</b>	<b>11</b>
<b>8.1.3</b>	<b>Measurement cell .....</b>	<b>11</b>
<b>8.1.4</b>	<b>Sample inspection, preparation, dilution, and concentration .....</b>	<b>12</b>
<b>8.2</b>	<b>Verification .....</b>	<b>13</b>
<b>8.2.1</b>	<b>Reference materials .....</b>	<b>13</b>
<b>8.2.2</b>	<b>Repeatability .....</b>	<b>13</b>
<b>8.2.3</b>	<b>Intermediate precision .....</b>	<b>13</b>
<b>8.2.4</b>	<b>Trueness .....</b>	<b>13</b>
<b>8.3</b>	<b>Sources of measurement error .....</b>	<b>14</b>
<b>8.3.1</b>	<b>Contamination of the current sample by the previous sample .....</b>	<b>14</b>
<b>8.3.2</b>	<b>Inappropriate sample preparation procedure .....</b>	<b>14</b>
<b>8.3.3</b>	<b>Inappropriate sample .....</b>	<b>14</b>
<b>8.3.4</b>	<b>Inappropriate liquid medium .....</b>	<b>14</b>
<b>8.3.5</b>	<b>Poor temperature stabilization .....</b>	<b>14</b>
<b>8.3.6</b>	<b>Condensation on the illuminated surfaces .....</b>	<b>14</b>
<b>8.3.7</b>	<b>Particles, fingerprints or scratches on the optical surfaces .....</b>	<b>14</b>
<b>8.3.8</b>	<b>Too large a potential applied .....</b>	<b>15</b>
<b>8.3.9</b>	<b>Incorrect entry of parameters by the operator .....</b>	<b>15</b>
<b>8.3.10</b>	<b>Air bubbles .....</b>	<b>15</b>

8.3.11	Cell coating damage .....	15
8.3.12	Calculating zeta-potential .....	15
8.3.13	Sample stability consideration .....	15
8.4	Test report .....	15
Annex A (informative) Electroosmosis within capillary cells .....		17
Bibliography .....		20