

# DIN EN ISO 2615:2024-10 (E)

## Analysis of natural gas - Biomethane - Determination of the content of compressor oil (ISO 2615:2024)

---

Contents	Page
European foreword.....	3
Foreword.....	4
Introduction.....	5
<b>1 Scope.....</b>	<b>6</b>
<b>2 Normative references.....</b>	<b>6</b>
<b>3 Terms and definitions.....</b>	<b>6</b>
<b>4 Principle.....</b>	<b>7</b>
<b>5 Chemicals and materials.....</b>	<b>7</b>
5.1 Compressor oils.....	7
5.2 Solvent.....	8
5.3 Sampling filters.....	8
5.4 Backup filters.....	8
5.5 Calibration standards.....	8
5.6 Nitrogen pure gas.....	8
5.7 Compressed liquid propane.....	8
<b>6 Apparatus.....</b>	<b>8</b>
6.1 Gas chromatograph.....	8
6.2 Capillary column.....	8
6.3 Ultrasonic bath.....	8
6.4 Pressurised fluid extraction apparatus.....	8
6.5 Filter housing.....	9
6.6 Rotary evaporator, concentration apparatus.....	9
6.7 Laboratory glassware.....	9
<b>7 Sampling.....</b>	<b>9</b>
7.1 Sampling apparatus.....	9
7.2 Sampling steps.....	9
<b>8 Methods to recover the oil from the buffer tank.....</b>	<b>10</b>
<b>9 Extraction procedure for coalescing filters.....</b>	<b>10</b>
9.1 Procedure 1: Ultrasonic extraction – nitrogen flush.....	10
9.2 Procedure 2: Pressurised fluid extraction.....	10
<b>10 Analysis.....</b>	<b>11</b>
10.1 GC/MS analysis.....	11
10.2 GC/FID analysis.....	11
<b>11 Calculations.....</b>	<b>11</b>
<b>12 Test report.....</b>	<b>12</b>
<b>Annex A (informative) Example of a GC/MS chromatogram obtained for a compressor oil (Figure A.1).....</b>	<b>13</b>
<b>Annex B (informative) Example of a GC/FID chromatogram obtained for a compressor oil (Figure B.1).....</b>	<b>14</b>
<b>Annex C (informative) Different types of oils.....</b>	<b>16</b>
<b>Bibliography.....</b>	<b>18</b>