

# ISO 5361:2016-09 (E)

## Anaesthetic and respiratory equipment - Tracheal tubes and connectors

---

<b>Contents</b>		<b>Page</b>
Foreword .....		iv
Introduction .....		v
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Normative references .....</b>	<b>1</b>
<b>3</b>	<b>Terms and definitions .....</b>	<b>2</b>
<b>4</b>	<b>*General requirements for tracheal tubes and tracheal tube connectors .....</b>	<b>5</b>
<b>4.1</b>	<b>Risk assessment .....</b>	<b>5</b>
<b>4.2</b>	<b>Safety .....</b>	<b>6</b>
<b>5</b>	<b>Specific requirements for tracheal tubes and tracheal tube connectors .....</b>	<b>6</b>
<b>5.1</b>	<b>Size designation .....</b>	<b>6</b>
<b>5.2</b>	<b>Dimensions .....</b>	<b>6</b>
<b>5.2.1</b>	<b>Tracheal tubes .....</b>	<b>6</b>
<b>5.2.2</b>	<b>Tracheal tube connectors .....</b>	<b>10</b>
<b>5.3</b>	<b>*Materials .....</b>	<b>13</b>
<b>5.4</b>	<b>Tracheal tube bevel .....</b>	<b>13</b>
<b>5.5</b>	<b>*Tracheal tube cuffs .....</b>	<b>13</b>
<b>5.6</b>	<b>Inflating system for cuffs .....</b>	<b>14</b>
<b>5.7</b>	<b>Curvature of the tube .....</b>	<b>14</b>
<b>5.8</b>	<b>*Radiopaque marker .....</b>	<b>16</b>
<b>5.9</b>	<b>*Kink resistance .....</b>	<b>16</b>
<b>6</b>	<b>Additional requirement for tracheal tubes with a Murphy eye .....</b>	<b>17</b>
<b>6.1</b>	<b>Size of the Murphy eye .....</b>	<b>17</b>
<b>6.2</b>	<b>Location of the Murphy eye .....</b>	<b>17</b>
<b>6.3</b>	<b>Resistance to kinking or collapse of the Murphy eye .....</b>	<b>17</b>
<b>6.4</b>	<b>Surface characteristics of the Murphy eye .....</b>	<b>17</b>
<b>7</b>	<b>Requirements for tracheal tubes with tracheal tube connectors supplied sterile .....</b>	<b>17</b>
<b>7.1</b>	<b>Sterility assurance .....</b>	<b>17</b>
<b>7.2</b>	<b>Packaging for tracheal tubes and tracheal tube connectors supplied sterile .....</b>	<b>18</b>
<b>8</b>	<b>Marking .....</b>	<b>18</b>
<b>8.1</b>	<b>Use of symbols .....</b>	<b>18</b>
<b>8.2</b>	<b>Marking on the tracheal tube .....</b>	<b>18</b>
<b>8.3</b>	<b>Marking on the tracheal tube individual pack or any insert .....</b>	<b>19</b>
<b>8.4</b>	<b>Marking on tracheal tube connectors .....</b>	<b>21</b>
<b>Annex A (informative)</b>	<b>Rationale .....</b>	<b>22</b>
<b>Annex B (normative)</b>	<b>Determination of cuff diameter .....</b>	<b>26</b>
<b>Annex C (normative)</b>	<b>Test method for cuffed tube collapse .....</b>	<b>27</b>
<b>Annex D (normative)</b>	<b>*Test method for cuff herniation .....</b>	<b>30</b>
<b>Annex E (informative)</b>	<b>Guidance on design of tracheal tube connectors .....</b>	<b>32</b>

<b>Annex F (informative) Hazard identification for risk assessment .....</b>	<b>33</b>
<b>Annex G (normative) *Test method for tracheal seal .....</b>	<b>36</b>
<b>Annex H (normative) Test method to determine kink resistance .....</b>	<b>39</b>
<b>Bibliography .....</b>	<b>41</b>