

DIN 13277:2022-05 (E)

Refrigerators and freezers for laboratory and medical applications - Terminology, requirements, testing

Contents	Page
Foreword	4
1 Scope	5
2 Normative references	6
3 Terms and definitions.....	6
4 Requirements	9
4.1 General	9
4.1.1 Gross volume	9
4.1.2 Usable volume	9
4.1.3 Installation dimensions	9
4.1.4 Interior space and interior fittings.....	9
4.1.5 Latching system for appliance door/flap and lid.....	9
4.1.6 Appliance door and lid seal.....	10
4.1.7 Installation and set-up.....	10
4.2 Safety requirements.....	10
4.2.1 Stability.....	10
4.2.2 Switch-on indicator	10
4.2.3 Usable-space temperature and product temperature	10
4.2.4 Monitoring of the product temperature	11
4.2.5 Temperature rise inside in case of power failure.....	13
4.2.6 Noise emission	13
4.2.7 Automatic defrosting	13
4.2.8 Appliance insulation / condensation formation	13
4.2.9 Relative duty cycle / needed capacity	13
4.3 Hygienic requirements.....	13
5 Testing	14
5.1 General	14
5.1.1 General test conditions	14
5.1.2 Dimensions.....	16
5.1.3 Gross volume	16
5.1.4 Usable volume	16
5.1.5 Loading capacity of the storage surfaces.....	16
5.1.6 Latching system for appliance door/flap and lid.....	16
5.1.7 Appliance door seal and lid seal.....	16
5.2 Safety requirements.....	17
5.2.1 General	17
5.2.2 Usable-space temperature and product temperature	17
5.2.3 Appliance insulation / condensation formation	24
5.2.4 Relative duty cycle / needed capacity.....	24
5.2.5 Standard energy consumption.....	25
6 Test report.....	25
7 Labelling.....	26

Annex A (informative) Exemplary method for product temperature testing of loaded blood bank refrigerators and loaded plasma storage appliances	27
A.1 Loaded appliance	27
A.1.1 General	27
A.1.2 Test bags	27
A.1.3 Loading of plasma storage appliances	27
A.1.4 Loading of blood bank refrigerators	29
A.1.5 Reference body	31
A.2 Testing the product temperature with loading	32
A.2.1 Testing the temperature distribution during a defrosting process for forced-air appliances	32
A.2.2 Testing the temperature distribution while opening the appliance door/lid/flap	32
A.2.3 Testing the temperature curve in case of power failure	32
Annex B (informative) Schematic overview of important terms for temperature testing	33
Bibliography	34

Figures

Figure 1 — Example of a suitable reference body set-up	12
Figure 2 — Partitions to reduce air circulation (top view)	15
Figure 3 — Overview of measuring planes	17
Figure 4 — Arrangement and position of measurement points	18
Figure 5 — Distances from storage surfaces and wall limits	19
Figure 6 — Test evaluation	20
Figure A.1 — Arrangement of test bags for loading upright plasma storage appliances	28
Figure A.2 — Arrangement of test bags for loading horizontal plasma storage appliances with and without appliance step	29
Figure A.3 — Arrangement of test bags for loading upright blood bank refrigerators	30
Figure A.4 — Arrangement of test bags for loading horizontal blood bank refrigerators with and without appliance step	31
Figure B.1 — Schematic overview of important terms for temperature testing	33

Tables

Table 1 — Storage applications of refrigerating and freezing storage appliances	5
Table 2 — Usable-space temperature range (measured with measuring body defined in 3.9) and product temperature range (measured with reference body defined in 3.16) of the refrigerating appliance in steady-state operation	11
Table 3 — List of all required tests and test requirements and related subclauses in this document	14
Table 4 — Limit temperature (measured in the reference body defined in 3.16) of the refrigerating appliance in case of power failure	23