

DIN EN 16159:2012-04 (E)

Animal feeding stuffs - Determination of selenium by hydride generation atomic absorption spectrometry (HGAAS) after microwave digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide)

Contents	Page
Foreword	3
1 Scope	4
2 Normative references	4
3 Principle	4
4 Reagents	5
5 Apparatus and equipment	6
6 Procedure	7
6.1 General	7
6.2 Preparation of the test solution	7
6.3 Measurement of the test solution	8
6.3.1 Pre-dilution of the test solution	8
6.3.2 Pre-reduction of the (pre-diluted) test solution	8
6.3.3 Spectrometer settings of the flow-injection hydride generation atomic absorption spectrometer (HGAAS)	9
6.3.4 HGAAS determination	9
7 Calculation	10
8 Precision	10
8.1 Introduction	10
8.2 General	11
8.3 Repeatability	11
8.4 Reproducibility	11
9 Test report	11
Annex A (informative) Results of the inter-laboratory test	12
Annex B (informative) Flowchart - Determination of selenium by hydride generation atomic absorption spectrometry (HGAAS) after microwave digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide)	13
Annex C (informative) Alternative digestion procedure with the same digestion efficiency: Acid digestion with a mixture of 65 % nitric acid and 70 % perchloric acid (7:3 by volume) at atmospheric pressure	14
C.1 Preparation of selenium calibration solutions	14
C.2 Preparation of the test solution	14
C.3 Measurement of the test solution	15
C.4 Ten rules for automated wet ashing with perchloric acid	15
Bibliography	16