

DIN EN 16274:2012-12 (E)

Methods for analysis of allergens - Quantification of suspected fragrance allergens in consumer products - Step 1: GC analysis of ready-to-inject sample

Contents		Page
Foreword		3
Introduction		4
1	Scope	5
2	Principle	5
3	Reagents	5
3.1	Solvents	5
3.2	Fragrance (suspected allergen) standards	5
3.3	Internal standards (ISTD)	7
4	Apparatus	7
4.1	Analytical Balance	7
4.2	GC-FID (for solvent or standard purity only)	7
4.3	GC-MS	7
4.4	GC capillary columns	7
5	Standard preparation and preservation	8
5.1	General	8
5.2	Standard preparation	8
5.2.1	General	8
5.2.2	Stock solution of all allergens (5 g/l)	8
5.2.3	Separate stock solutions (carbonyl / non carbonyl compounds) (10 g/l)	9
5.2.4	Internal standard solution	9
5.2.5	Working solutions	9
5.2.6	Calibration solution	9
6	Procedure	10
6.1	General	10
6.2	Chromatographic conditions	10
6.3	MS conditions	10
6.3.1	General	10
6.3.2	SCAN mode	10
6.3.3	SIM mode	11
6.4	Calibration	13
7	Sample analysis	13
8	Data treatment and calculation of results	14
8.1	Identification of allergens	14
8.2	Quantification of allergens	14
8.3	Assessment of the analytical measurement	15
8.3.1	General	15
8.3.2	Examination of the Q-values	15
8.3.3	Maximum permitted tolerances	15
9	Test report	16
Annex A (informative)	Column performances	17

Annex B (informative) SIM windows	19
Annex C (informative) Example of chromatographic separation	20
Annex D (normative) Decisional tree	22
Bibliography	23