

ISO 27971:2015-06 (E)

Cereals and cereal products - Common wheat (*Triticum aestivum* L.) - Determination of alveograph properties of dough at constant hydration from commercial or test flours and test milling methodology

Contents		Page
Foreword		v
Introduction		vi
1	Scope	1
2	Normative references	1
3	Principle	1
4	Reagents	1
5	Apparatus	2
6	Sampling	9
7	Preparation of the wheat for laboratory milling	10
7.1	Cleaning the laboratory sample	10
7.2	Test portion	10
7.3	Wheat moisture content determination	10
7.4	Wheat preparation	10
7.4.1	General	10
7.4.2	Wheat with initial moisture content between 13 % and 15 % (one-stage moistening)	10
7.4.3	Wheat with a moisture content less than 13 % (two-stage moistening)	10
7.4.4	Wheat with a moisture content greater than 15 % (preliminary drying followed by moistening, as described above)	11
8	Laboratory milling	11
8.1	General	11
8.2	Milling procedure	11
8.2.1	Breaking	11
8.2.2	Reduction	11
8.2.3	Flour homogenization	12
8.2.4	Storage of the flour	12
8.3	Expression of milling results	12
9	Preparation and alveograph test	13
9.1	Preliminary checks	13
9.2	Preliminary operations	14
9.3	Kneading	15
9.4	Preparation of dough test pieces	16
9.5	Alveograph test	20
9.5.1	Initial preparation	20
9.5.2	First operation: Adjusting the dough test piece	21
9.5.3	Second operation: biaxial extension	22
9.6	Expression of alveograph test results	23
9.6.1	General	23
9.6.2	Maximum pressure parameter, P	23
9.6.3	Mean abscissa at rupture, L	23

9.6.4	Swelling index, G	24
9.6.5	Elasticity index	24
9.6.6	Curve configuration ratio, P/L	24
9.6.7	Deformation work, W	24
10	Precision	24
10.1	Interlaboratory tests	24
10.2	Repeatability limits	25
10.2.1	Commercial flour: limits established by the interlaboratory test	25
10.2.2	Flour obtained from laboratory milling	25
10.3	Reproducibility limits	25
10.3.1	Commercial flour: Limits established by the proficiency tests	25
10.3.2	Flour obtained from laboratory milling	26
10.4	Uncertainty	26
11	Test report	26
Annex A (informative)	Characteristics of the Chopin-Dubois CD1 mill	27
Annex B (normative)	Quantity of water to be added to wheat for conditioning	29
Annex C (informative)	Sample milling sheet	31
Annex D (informative)	Conversion table from L to G	32
Annex E (informative)	Interlaboratory and proficiency test data for commercial flours	34
Annex F (informative)	Interlaboratory data for laboratory milled flour	44
Annex G (informative)	Routine maintenance instructions for the alveograph	61
Annex H (informative)	Assessment of proteolytic activity in wheat (<i>T. aestivum</i> L.) or flour	63
Bibliography	65