

# DIN EN 16877:2017-03 (E)

Animal feeding stuffs - Methods of sampling and analysis - Determination of T-2 and HT-2 toxins, Deoxynivalenol and Zearalenone, in feed materials and compound feed by LC-MS

---

<b>Contents</b>		<b>Page</b>
European foreword .....		4
Introduction .....		5
1	Scope .....	6
2	Normative references .....	6
3	Principle .....	6
4	Reagents .....	6
5	Apparatus .....	8
6	Procedures .....	9
6.1	Sample preparation .....	9
6.2	Extraction .....	9
6.3	Test solution .....	10
6.4	Spiking procedure .....	10
7	Measurements .....	11
7.1	General .....	11
7.2	LC conditions .....	11
7.3	MS conditions .....	11
7.4	Batch composition .....	11
7.5	Peak identification .....	11
7.6	Determination of DON, HT2, T2, and ZON in calibration or and test solutions .....	11
7.7	Calibration .....	11
8	Determination of mass fraction .....	12
9	Precision .....	13
9.1	Interlaboratory study .....	13
9.2	Repeatability .....	13
9.2.1	General .....	13
9.2.2	HT-2 toxin .....	13
9.2.3	T-2 toxin .....	13
9.2.4	DON .....	13
9.2.5	ZON .....	13
9.3	Reproducibility .....	13
9.3.1	General .....	13
9.3.2	HT-2 toxin .....	14
9.3.3	T-2 toxin .....	14
9.3.4	DON .....	14
9.3.5	ZON .....	14
10	Test report .....	14
Annex A (informative) Precision data .....		15

<b>Annex B (informative) Examples .....</b>	<b>20</b>
<b>B.1 Example 1 .....</b>	<b>20</b>
<b>B.1.1 General .....</b>	<b>20</b>
<b>B.1.2 LC conditions .....</b>	<b>20</b>
<b>B.1.3 MS conditions .....</b>	<b>21</b>
<b>B.2 Example 2 .....</b>	<b>21</b>
<b>B.2.1 General .....</b>	<b>21</b>
<b>B.2.2 LC conditions .....</b>	<b>22</b>
<b>B.2.3 MS conditions .....</b>	<b>22</b>
<b>B.3 Example 3 .....</b>	<b>23</b>
<b>B.3.1 General .....</b>	<b>23</b>
<b>B.3.2 LC conditions .....</b>	<b>24</b>
<b>B.3.3 MS conditions .....</b>	<b>24</b>
<b>Annex C (informative) Examples of chromatograms according to the settings of the examples in Annex B .....</b>	<b>26</b>
<b>Bibliography .....</b>	<b>31</b>