

# DIN EN 17212:2019-11 (Englisch)

## Animal Feeding stuffs - Methods of sampling and analysis - Determination of melamine and cyanuric acid content by liquid chromatographic method with mass spectrometric detection (LC-MS/MS)

---

<b>Contents</b>		<b>Page</b>
European foreword .....		4
Introduction .....		5
1	Scope .....	6
2	Normative references .....	6
3	Terms and definitions .....	7
4	Principle .....	11
5	Reagents .....	11
5.1	General .....	11
5.2	Chemicals .....	12
5.3	Standard solutions .....	13
5.3.6	Cyanuric acid stock solution 1, c = 1 mg/ml .....	14
6	Apparatus .....	15
7	Preparation of sample .....	16
7.1	General .....	16
7.2	Laboratory sample .....	16
7.3	Test material .....	16
7.4	Sample material .....	16
8	Procedure .....	17
8.1	General .....	17
8.2	Screening method .....	17
8.2.1	Sample extraction .....	17
8.2.2	Sample dilution .....	17
8.2.3	Injection solutions for LC-MS/MS .....	17
8.3	Confirmatory method .....	17
8.3.1	Addition of internal standard .....	17
8.3.2	Sample extraction .....	18
8.3.3	Sample dilution and standard addition .....	18
8.3.4	Injection solution for LC-MS/MS .....	19
8.4	Quality control measures for screening method .....	20
8.4.1	Blank extract .....	20
8.4.2	Control standards .....	20
8.5	Quality control measures for confirmatory method .....	20
8.5.1	Blank extract .....	20
8.5.2	Control matrix sample .....	20
8.5.3	Control standards .....	21
9	HPLC-MS/MS analysis .....	21
9.1	General .....	21
9.2	HPLC operating conditions .....	21
9.3	Determination of melamine and cyanuric acid in sample test solutions .....	21

9.3.1	Example of an analysis sequence for screening method .....	21
9.3.2	Sample tested positive in the screening method .....	22
9.3.3	Example of an analysis sequence for confirmatory method .....	22
10	Calculation .....	22
10.1	Peak identification .....	22
10.2	Quantification (confirmatory method) .....	22
11	Precision .....	24
11.1	Collaborative trial .....	24
11.2	Repeatability limit .....	24
11.3	Reproducibility limit .....	24
12	Test report .....	25
Annex A (informative) Examples for suitable HPLC-MS/MS conditions .....		26
A.1	Example of suitable HPLC conditions .....	26
A.2	Example suitable for SCIEX API 4000 Q-Trap or API 4000 .....	27
Annex B (informative) Results of the international collaborative trial .....		29
Annex C (informative) Results of the in-house validation of the method .....		32
Annex D (informative) Typical chromatograms .....		35
Bibliography .....		37