

DIN EN 16274:2021-11 (E)

Method for analysis of allergens - Quantification of an extended list of 57 suspected allergens in ready to inject fragrance materials by gas chromatography mass spectrometry

Contents		Page
European foreword		4
Introduction		5
1	Scope	6
2	Normative references	6
3	Terms and definitions	6
4	Principle	6
5	Reagents	7
5.1	Solvents	7
5.1.1	Methyl pivalate	7
5.1.2	Tert-Butyl Methyl Ether (MTBE)	7
5.2	Reference samples (suspected allergens)	7
5.2.1	General	7
6	Apparatus	11
6.1	Gas chromatograph equipped with flame ionization detector (GC-FID)	11
6.2	Gas chromatograph coupled to a mass spectrometer (GC-MS)	11
6.2.1	General	11
6.2.2	GC-MS system	11
6.3	Capillary columns for GC	11
6.4	Analytical balance	12
7	Stock and sample solutions - preparation and storage	12
7.1	General information	12
7.1.1	Choice of the solvent	12
7.1.2	Miscellaneous	13
7.2	Preparation of stock solutions from reference samples	13
7.3	Preparation of internal standard solutions and calibration solutions	13
7.3.1	General	13
7.3.2	Stock Solutions of Internal Standards (1g/kg equivalent)	13
7.3.3	Calibration Solutions Concentrations	13
7.3.4	Preparation of Calibration solutions	14
8	Mass spectrometer acquisition conditions	16
8.1	Establishment of retention times and SIM chromatograph windows	16
8.2	SIM Window - Set up and criteria	17
8.3	GC-MS scan mode verification	17
9	Sample analysis	17
9.1	General	17
9.2	Sample preparation for analysis	17
9.3	Sequence	18
9.3.1	Blank	18
9.3.2	Calibration solutions	18

9.3.3	Check solution	18
9.3.4	Samples	19
10	Data validation and treatment	19
10.1	General	19
10.2	Examination of Q values	19
10.3	Relative intensity of characteristic ions in SIM	20
10.4	Relative ion intensity using scan data	20
10.5	Data verification scheme and reporting of final concentration	21
10.6	GC-MS scan mode verification	21
10.7	Mass spectrum evaluation in SCAN mode	21
11	Test report	21
Annex A (informative) Determination of calibrant and reference sample purity		22
Annex B (informative) GC capillary column parameters		29
Annex C (informative) SIM Ions for SIM or SIM-SCAN Use		30
Annex D (informative) Example of SIM windows		37
Annex E (informative) Example of chromatograms		46
Annex F (normative) Decisional tree for quantification of suspected allergens		48
Annex G (informative) Preparation of stock solutions from reference samples		49
Annex H (informative) Calibration methods and approach		54
Annex I (informative) Quantification of allergens		56
Annex J (informative) Suspected allergens selected as analytes targets - Rationale		59
Annex K (informative) Other general information		69
Bibliography		70