

# ISO 28218:2010-10 (E)

## Radiation protection - Performance criteria for radiobioassay

---

<b>Contents</b>		<b>Page</b>
Foreword .....		iv
Introduction .....		v
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Normative references .....</b>	<b>2</b>
<b>3</b>	<b>Terms and definitions .....</b>	<b>2</b>
<b>4</b>	<b>Symbols .....</b>	<b>6</b>
<b>5</b>	<b>Performance measures .....</b>	<b>7</b>
<b>5.1</b>	<b>Decision threshold (<math>y^*</math>) and detection limit (<math>y\#</math>) .....</b>	<b>7</b>
<b>5.2</b>	<b>Relative bias and bias performance criteria .....</b>	<b>10</b>
<b>5.3</b>	<b>Repeatability performance criteria .....</b>	<b>12</b>
<b>6</b>	<b>Performance criteria for in vivo radiobioassay .....</b>	<b>12</b>
<b>6.1</b>	<b>General .....</b>	<b>12</b>
<b>6.2</b>	<b>Responsibilities of the customer that could impact the service laboratory's performance .....</b>	<b>12</b>
<b>6.3</b>	<b>Service laboratory criteria .....</b>	<b>12</b>
<b>6.4</b>	<b>Identification of radionuclides .....</b>	<b>13</b>
<b>6.5</b>	<b>Quantification .....</b>	<b>13</b>
<b>6.6</b>	<b>Reporting results .....</b>	<b>14</b>
<b>6.7</b>	<b>Records retention .....</b>	<b>15</b>
<b>7</b>	<b>Performance criteria for in vitro radiobioassay .....</b>	<b>15</b>
<b>7.1</b>	<b>General .....</b>	<b>15</b>
<b>7.2</b>	<b>Responsibilities of the customer that could impact the service laboratory's performance .....</b>	<b>15</b>
<b>7.3</b>	<b>Analytical methodology .....</b>	<b>16</b>
<b>7.4</b>	<b>Reporting results .....</b>	<b>17</b>
<b>7.5</b>	<b>Records retention .....</b>	<b>17</b>
<b>8</b>	<b>Quality assurance and quality control for radiobioassay laboratories .....</b>	<b>18</b>
<b>8.1</b>	<b>General .....</b>	<b>18</b>
<b>8.2</b>	<b>Quality assurance .....</b>	<b>18</b>
<b>8.3</b>	<b>Quality assurance plan .....</b>	<b>19</b>
<b>8.4</b>	<b>Quality control .....</b>	<b>20</b>
<b>9</b>	<b>Performance testing programme .....</b>	<b>21</b>
<b>9.1</b>	<b>General .....</b>	<b>21</b>
<b>9.2</b>	<b>In vivo radiobioassay .....</b>	<b>22</b>
<b>9.3</b>	<b>In vitro radiobioassay .....</b>	<b>25</b>
	<b>Annex A (informative) Detection limit -- Models for applications .....</b>	<b>29</b>
	<b>Annex B (informative) Detection limit -- Application examples .....</b>	<b>32</b>
	<b>Bibliography .....</b>	<b>45</b>