

# DIN 38407-42:2011-03 (E)

German standard methods for the examination of water, waste water and sludge -  
Jointly determinable substances (group F) - Part 42: Determination of selected  
polyfluorinated compounds (PFC) in water - Method using high performance liquid  
chromatography and mass spectrometric detection (HPLC/MS-MS) after solid-liquid  
extraction (F 42)

---

<b>Contents</b>	<b>Page</b>
Introduction .....	5
1 Scope .....	6
2 Normative references .....	7
3 Terms and definitions .....	7
4 Principle .....	7
5 Interferences .....	8
5.1 General .....	8
5.2 Interferences encountered during extraction and processing of extracts .....	8
5.3 Interferences encountered during high performance liquid chromatography and mass spectrometry .....	8
6 Designation .....	9
7 Reagents .....	9
8 Apparatus .....	11
9 Sampling .....	12
10 Procedure .....	12
10.1 General .....	12
10.2 Sample preparation .....	12
10.3 Extraction .....	12
10.4 High performance liquid chromatography (HPLC) .....	13
10.5 Detection .....	14
10.6 Blank value measurements .....	14
11 Calibration .....	15
11.1 Principles .....	15
11.2 Calibration using an external standard .....	16
11.3 Calibration using an internal standard .....	17
11.4 Calibration check .....	18
12 Determination of recoveries .....	19
12.1 Recoveries of the method .....	19
12.2 Recovery rates of internal standards .....	20
13 Evaluation .....	21
13.1 General .....	21
13.2 Verification of individual substances .....	21
13.3 Calculation of the individual result .....	22

14	Expression of results .....	23
15	Test report .....	23
16	Performance data .....	23
Annex A (informative) Examples of sorbents, working conditions and recovery rates .....		27
A.1	Examples of sorbents and working conditions suitable for solid-phase extraction .....	27
Annex B (informative) Examples of suitable HPLC columns and chromatograms .....		34
B.1	Chromatographic conditions for the chromatogram shown in Figure B.1 .....	34
B.2	Chromatographic conditions for the chromatogram shown in Figure B.2 .....	35
B.3	Chromatographic conditions for the chromatogram shown in Figure B.3 .....	36
B.4	Chromatographic conditions for the chromatogram shown in Figure B.4 .....	37
B.5	Chromatographic conditions for the chromatogram shown in Figure B.5 .....	38
Annex C (informative) Examples of selected diagnostic ions for identification and quantification ...		39
Annex D (informative) Examples for the expansion of the method .....		41
Annex E (informative) Explanatory notes .....		42
Bibliography .....		43
Figures Figure B.1 -- Chromatographic separation, example 1 .....		34
Figure B.2 -- Chromatographic separation, example 2 .....		35
Figure B.3 -- Chromatographic separation, example 3 .....		36
Figure B.4 -- Chromatographic separation, example 4 .....		37
Figure B.5 -- Chromatographic separation, example 5 .....		38
Figure C.1 -- Example of a MS chromatogram of a surface water sample (extract shows branched and unbranched perfluorinated carboxylic acids) .....		40
Figure C.2 -- Example of a MS chromatogram of a surface water sample (extract shows branched and unbranched perfluorinated sulfonic acids) .....		40
Tables Table 1 -- Substances whose determination has been tested using this method .....		6
Table 2 -- Meaning of the indices .....		16
Table 3 -- Example of the assignment of internal standard substances to the analytes .....		17
Table 4 -- Performance data for HPLC/MS-MS measurement .....		24
Table 5 -- Performance data for drinking water .....		25
Table 6 -- Performance data for ground water .....		25
Table 7 -- Performance data for surface water .....		26
Table 8 -- Performance data for treated waste water .....		26
Table A.1 -- Examples of recovery rates - Ultra pure water (as in 7.2) .....		30
Table A.2 -- Examples of recovery rates - Drinking water .....		31

<b>Table A.3 -- Examples of recovery rates - Surface water .....</b>	<b>32</b>
<b>Table A.4 -- Examples of recovery rates - Treated waste water (effluents) .....</b>	<b>33</b>
<b>Table C.1 -- Selected diagnostic ions (target compounds) .....</b>	<b>39</b>
<b>Table C.2 -- Selected diagnostic ions (internal standards) .....</b>	<b>39</b>