

DIN 38407-42:2011-03 (E)

German standard methods for the examination of water, waste water and sludge -
Jointly determinable substances (group F) - Part 42: Determination of selected
polyfluorinated compounds (PFC) in water - Method using high performance liquid
chromatography and mass spectrometric detection (HPLC/MS-MS) after solid-liquid
extraction (F 42)

| Contents | Page |
|---|-------------|
| Introduction | 5 |
| 1 Scope | 6 |
| 2 Normative references | 7 |
| 3 Terms and definitions | 7 |
| 4 Principle | 7 |
| 5 Interferences | 8 |
| 5.1 General | 8 |
| 5.2 Interferences encountered during extraction and processing of extracts | 8 |
| 5.3 Interferences encountered during high performance liquid chromatography and mass spectrometry | 8 |
| 6 Designation | 9 |
| 7 Reagents | 9 |
| 8 Apparatus | 11 |
| 9 Sampling | 12 |
| 10 Procedure | 12 |
| 10.1 General | 12 |
| 10.2 Sample preparation | 12 |
| 10.3 Extraction | 12 |
| 10.4 High performance liquid chromatography (HPLC) | 13 |
| 10.5 Detection | 14 |
| 10.6 Blank value measurements | 14 |
| 11 Calibration | 15 |
| 11.1 Principles | 15 |
| 11.2 Calibration using an external standard | 16 |
| 11.3 Calibration using an internal standard | 17 |
| 11.4 Calibration check | 18 |
| 12 Determination of recoveries | 19 |
| 12.1 Recoveries of the method | 19 |
| 12.2 Recovery rates of internal standards | 20 |
| 13 Evaluation | 21 |
| 13.1 General | 21 |
| 13.2 Verification of individual substances | 21 |
| 13.3 Calculation of the individual result | 22 |

| | | |
|--|---|----|
| 14 | Expression of results | 23 |
| 15 | Test report | 23 |
| 16 | Performance data | 23 |
| Annex A (informative) Examples of sorbents, working conditions and recovery rates | | 27 |
| A.1 | Examples of sorbents and working conditions suitable for solid-phase extraction | 27 |
| Annex B (informative) Examples of suitable HPLC columns and chromatograms | | 34 |
| B.1 | Chromatographic conditions for the chromatogram shown in Figure B.1 | 34 |
| B.2 | Chromatographic conditions for the chromatogram shown in Figure B.2 | 35 |
| B.3 | Chromatographic conditions for the chromatogram shown in Figure B.3 | 36 |
| B.4 | Chromatographic conditions for the chromatogram shown in Figure B.4 | 37 |
| B.5 | Chromatographic conditions for the chromatogram shown in Figure B.5 | 38 |
| Annex C (informative) Examples of selected diagnostic ions for identification and quantification ... | | 39 |
| Annex D (informative) Examples for the expansion of the method | | 41 |
| Annex E (informative) Explanatory notes | | 42 |
| Bibliography | | 43 |
| Figures Figure B.1 -- Chromatographic separation, example 1 | | 34 |
| Figure B.2 -- Chromatographic separation, example 2 | | 35 |
| Figure B.3 -- Chromatographic separation, example 3 | | 36 |
| Figure B.4 -- Chromatographic separation, example 4 | | 37 |
| Figure B.5 -- Chromatographic separation, example 5 | | 38 |
| Figure C.1 -- Example of a MS chromatogram of a surface water sample (extract shows branched and unbranched perfluorinated carboxylic acids) | | 40 |
| Figure C.2 -- Example of a MS chromatogram of a surface water sample (extract shows branched and unbranched perfluorinated sulfonic acids) | | 40 |
| Tables Table 1 -- Substances whose determination has been tested using this method | | 6 |
| Table 2 -- Meaning of the indices | | 16 |
| Table 3 -- Example of the assignment of internal standard substances to the analytes | | 17 |
| Table 4 -- Performance data for HPLC/MS-MS measurement | | 24 |
| Table 5 -- Performance data for drinking water | | 25 |
| Table 6 -- Performance data for ground water | | 25 |
| Table 7 -- Performance data for surface water | | 26 |
| Table 8 -- Performance data for treated waste water | | 26 |
| Table A.1 -- Examples of recovery rates - Ultra pure water (as in 7.2) | | 30 |
| Table A.2 -- Examples of recovery rates - Drinking water | | 31 |

| | |
|--|-----------|
| Table A.3 -- Examples of recovery rates - Surface water | 32 |
| Table A.4 -- Examples of recovery rates - Treated waste water (effluents) | 33 |
| Table C.1 -- Selected diagnostic ions (target compounds) | 39 |
| Table C.2 -- Selected diagnostic ions (internal standards) | 39 |