

DIN EN 13946:2014-07 (E)

Water quality - Guidance for the routine sampling and preparation of benthic diatoms from rivers and lakes

Contents		Page
Foreword		3
Introduction		4
1	Scope	5
2	Terms and definitions	5
3	Principle	6
4	Equipment	6
4.1	Equipment for field sampling	6
4.2	Laboratory equipment	6
5	Reagents	6
6	Procedure	8
6.1	Choice of substratum	8
6.2	Sample site selection	8
6.3	Sampling methods	8
6.3.1	Moveable natural hard surfaces	8
6.3.2	Method for sampling vertical man-made surfaces in situ	9
6.3.3	Use of artificial substrata	10
6.3.4	Sample collection from submerged macrophytes and macroalgae	10
6.3.5	Sample collection from emergent macrophytes	10
6.4	Preparation prior to microscopic examination	11
6.4.1	Preservation and preliminary laboratory treatment	11
6.4.2	Methods for cleaning diatoms	11
6.4.3	Preparation of permanent slides	11
Annex A (informative)	Methods for cleaning diatoms for microscopic examination	13
A.1	General	13
A.2	Method 1: Hot hydrogen peroxide	13
A.2.1	Apparatus	13
A.2.2	Reagents	13
A.2.3	Procedure	13
A.3	Method 2: Cold hydrogen peroxide	14
A.3.1	Apparatus and reagents	14
A.3.2	Procedure	14
A.4	Method 3: Hot hydrogen peroxide with potassium dichromate	14
A.4.1	Apparatus	14
A.4.2	Reagents	14
A.4.3	Procedure	14
A.5	Method 4: Cold acid (permanganate) method of cleaning	15
A.5.1	Apparatus	15
A.5.2	Reagents	15
A.5.3	Procedure	15
Bibliography		17