

ISO 19040-2:2018 (E)

Water quality — Determination of the estrogenic potential of water and waste water — Part 2: Yeast estrogen screen (A-YES, *Arxula adenivorans*)

Contents

	Foreword
1	Scope
2	Normative references
3	Terms and definitions
4	Principle
5	Interferences
6	Apparatus and materials
7	Reagents, media and test strains
8	Sampling and samples
8.1	General
8.2	Bottles and material for sampling
8.3	Bottles and material pre-cleaning
8.4	Sampling procedure
8.5	Transport of samples
8.6	Pretreatment of samples
8.7	Storage of samples
9	Procedure
9.1	Test set up
9.1.1	Preparation of the reference dilution series
9.1.2	Reactivation of the yeast
9.1.3	Negative control
9.1.4	Blank replicate
9.1.5	Sample dilution
9.1.6	Field blank
9.1.7	Plate set up
9.1.8	Inoculation of the test plate
9.1.8.1	Preparation of test sample specific inoculation medium
9.1.8.2	Inoculation of wells without test sample
9.1.8.3	Inoculation of wells with test sample
9.2	Measurement
9.2.1	Measurement of the reporter gene activity
9.2.2	Measurement of the cell density
9.3	Calculations
9.3.1	Background correction
9.3.2	Calculation of the relative growth
9.3.3	Calculations for assessment of sample blanks
9.3.3.1	Calculations for assessment of sample coloration and artificial cleavage of substrate
9.3.3.2	Calculations for assessment of sample turbidity
9.3.4	Calculation of the reporter gene induction
10	Validity criteria
11	Assessment criteria

12	Test report
13	Verification
Annex A (informative) Plate set up	
A.1	General
A.2	Scheme of plate set up
Annex B (informative) Lyophilization of <i>Arxula adenivorans</i> cells	
B.1	General
B.2	Reagents
B.3	Media
B.3.1	Glucose solution
B.3.2	Sorbitol solution
B.3.3	PBS buffer
B.3.4	Cryo medium
B.3.5	Lysin solution
B.3.6	Yeast minimal medium (YMM) with glucose
B.3.7	Yeast minimal medium (YMM) with sorbitol
B.4	Cultivation of <i>Arxula adenivorans</i> on solid agar
B.5	Cultivation of starter and main culture
B.6	Preparing of lyophilized <i>Arxula adenivorans</i> cells
Annex C (informative) Scheme of test principle	
Annex D (informative) Test set up for chemicals and extracts	
D.1	General
D.2	Extraction of water samples
D.3	Test with diluted organic solutions or extracts
Annex E (informative) Preparation of dilution series	
Annex F (informative) Performance data	
F.1	Validation study
F.1.1	Samples and participating laboratories
F.1.2	Results for lowest ineffective dilution (LID) and 17 β -estradiol equivalent concentrations (EEQ)
F.1.3	Trueness evaluation
F.2	Sensitivity for selected compounds
Annex G (informative) Statistical assessment	
Annex H (informative) Calculation of estradiol equivalents	
H.1	General
H.2	Modelling of dose response relationship
H.3	Calculation of estradiol equivalents for test samples
H.4	Reporting of estradiol equivalents for test samples
Annex I (informative) Alternative test design for EEQ determination	
Annex J (informative) Measurement of the lowest ineffective dilution (LID) of a waste water — A simplified evaluation for testing of waste water	
J.1	General
J.2	Principle
J.3	Preparation of dilutions for LID assessment
J.4	Test for LID
J.5	Assessment of results — LID value, effluents
J.6	Documentation of results
Annex K (informative) Example for statistical evaluation	
K.1	General
K.2	Calculations
K.2.1	Background correction

- K.2.2** Calculation of the relative growth
- K.2.3** Calculations for assessment of sample blanks
- K.2.3.1** Sample coloration
- K.2.3.2** Artificial cleavage of substrate
- K.2.4** Calculations for assessment of sample turbidity
- K.2.5** Calculation of the report gene induction
- K.2.5.1** Calculation of the corrected absorbance for sample dilutions
- K.2.5.2** Calculation of the corrected absorbance for dilutions of the reference compound and negative control
- K.2.6** Calculation of induction rates

Page count: 55