

# ISO 19040-3:2018 (E)

## Water quality — Determination of the estrogenic potential of water and waste water — Part 3: In vitro human cell-based reporter gene assay

---

### Contents

	Foreword
1	Scope
2	Normative references
3	Terms and definitions
4	Interferences
5	Principle
6	Apparatus and materials
7	Reagents, cells and media
8	Sampling and samples
8.1	General
8.2	Bottles and material for sampling
8.3	Bottles and material pre-cleaning
8.4	Sampling procedure
8.5	Transport of samples
8.6	Pretreatment of sample
8.7	Storage of samples
9	Procedure
9.1	Cell culture maintenance
9.1.1	Freezing cells
9.1.2	Starting a new cell culture
9.1.3	Culturing cells
9.2	Human cell reporter gene assay test procedure
9.2.1	Seeding the cells (day 1)
9.2.2	Preparation of the E2-reference (day 2)
9.2.3	Preparation of the sample dilutions (day 2)
9.2.4	Field blank
9.2.5	Exposing the cells (day 2)
9.2.6	Harvesting the cells (day 3)
9.2.7	Measurement of luminescence (day 3)
9.3	Data analysis
9.3.1	Calculation of the reporter gene induction
9.3.2	Calculation of the percentage of maximum response
9.3.3	Calculation of the dose-response curve
10	Validity criteria
10.1	Validity criteria for the assay
10.2	Validity criteria for samples
11	Assessment criteria
12	Test report
Annex A	(informative) Settings of the luminometer

**Annex B (informative) Plate setup**

**Annex C (informative) Bioassay characteristics and details**

**Annex D (informative) Test set up for chemicals and extracts**

- D.1 General**
- D.2 Extraction of water samples**
- D.3 Test with organic solutions or extracts**
- D.4 Data from literature**

**Annex E (informative) Preparation of dilution series**

**Annex F (informative) Performance data**

- F.1 Design of the interlaboratory trial**
  - F.1.1 General**
  - F.1.2 Description of samples and participating laboratories**
  - F.1.3 Summary of 17 $\beta$ -estradiol equivalent concentrations (EEQ)**
  - F.1.4 Summary of lowest ineffective dilutions LID**
  - F.1.5 Trueness of results**

**Annex G (informative) Statistical assessment**

**Annex H (informative) Calculation of 17 $\beta$ -estradiol equivalents**

- H.1 General**
- H.2 Modelling of dose response relationship**
- H.3 Calculation of 17 $\beta$ -estradiol equivalents for test samples**
- H.4 Reporting of 17 $\beta$ -estradiol equivalents for test samples**

**Annex I (informative) Measurement of the lowest ineffective dilution (LID) of waste water — A simplified evaluation for testing of waste water**

- I.1 General**
- I.2 Principle**
- I.3 Preparation of dilutions for LID assessment**
- I.4 Test for LID**
- I.5 Assessment of results — LID value. effluents**
- I.6 Reporting of results**

**Page count: 40**