

ISO 17858:2007-02 (E)

Water quality - Determination of dioxin-like polychlorinated biphenyls - Method using gas chromatography/mass spectrometry

Contents		Page
Foreword		vi
Introduction		vii
1	Scope	1
2	Normative references	1
3	Terms, definitions and abbreviated terms	2
3.1	Terms and definitions	2
3.2	Abbreviated terms	3
4	Principle	4
4.1	Spiking and extraction	4
4.2	Clean-up	4
4.3	Concentration	4
4.4	Identification	5
4.5	Quantification	5
4.6	Analytical quality	5
5	Contamination and interferences	5
6	Reagents and standards	6
7	Apparatus and materials	10
7.1	Sampling equipment for discrete sampling	10
7.2	Equipment for sample preparation	11
7.3	Extraction apparatus	11
7.4	Filtration apparatus	12
7.5	Clean-up apparatus	12
7.6	Concentration apparatus	13
7.7	Other equipment	13
8	Sample collection, preservation, storage and holding times	14
9	Quality assurance (QA)/quality control (QC)	14
9.1	General	14
9.2	Initial precision and recovery (IPR)	15
9.3	Spiking	15
9.4	Recovery of labelled compounds assessment	16
9.5	Method blanks	16
9.6	QC check sample	16
9.7	Method precision	16
10	Calibration	17
10.1	Operating conditions	17
10.2	Mass spectrometer (MS) resolution	17
10.3	Ion abundance ratios, minimum levels, signal-to-noise ratios, and absolute retention times	17
10.4	Retention time	18
10.5	Isomer specificity	18

10.6	Calibration by isotope dilution	18
10.7	Calibration by internal standard	19
10.8	Combined calibration	19
11	Sample preparation	20
11.1	General	20
11.2	Determination of percent suspended solids	20
11.3	Preparation of aqueous samples containing 1 % of suspended solids or less	21
12	Extraction and concentration	22
12.1	Separatory funnel extraction of filtrates and of aqueous samples that are visibly absent of particles	22
12.2	Solid-phase extraction (SPE) of samples containing less than 1 % suspended solids	22
12.3	Soxhlet extraction of filters and/or disks	23
12.4	Back-extraction with acid and base	24
12.5	Macro-concentration	24
12.6	Micro-concentration and solvent exchange	26
13	Extract clean-up	26
13.1	General	26
13.2	Gel permeation chromatography (GPC)	27
13.3	Silica clean-up	28
13.4	Alumina clean-up	28
13.5	Carbon column	29
13.6	High performance liquid chromatography (HPLC)	29
13.7	Florisil clean-up	30
13.8	Silver nitrate/silica column	31
14	HRGC/HRMS analysis	31
15	System and laboratory performance	31
15.1	General	31
15.2	MS resolution	31
15.3	Calibration verification	31
15.4	GC resolution	32
15.5	Blank	32
16	Qualitative determination	32
17	Quantitative determination	32
17.1	Isotope dilution quantification	32
17.2	Internal standard quantification and labelled-compound recovery	33
17.3	Concentration in sample	34
17.4	Results and reporting	35
17.5	Toxic equivalents (TEQ)	35
18	Analysis of complex samples	36
18.1	General	36
18.2	Recovery of labelled compounds	36
19	Pollution prevention	36
20	Waste management	37
21	Precision	37
	Annex A (informative) Example chromatograms	45
	Annex B (informative) Use of HRGC/LRMS	47
	Annex C (informative) Precision data	50

Bibliography	54
Table 1 -- Dioxin-like PCBs determined by this method	38
Table 2 -- Suggested quantification relationships	39
Table 3 -- Suggested calibration standard concentrations	40
Table 4 -- Suggested concentration of dioxin-like PCBs in stock and spiking solutions	41
Table 5 -- Typical GC columns and temperature programs	42
Table 6 -- Examples of toxic equivalent factors	43
Table 7 -- Congener function groups and ions	44
Table B.1 -- TetraCBs	49
Table B.2 -- PentaCBs	49
Table B.3 -- HexaCBs	49
Table B.4 -- HeptaCBs	49
Table C.1 -- Spiking amounts transferred to sample bottles	50
Table C.2 -- Samples 1 and 2 (fortified industrial effluent) -- Statistical summary	51
Table C.3 -- Sample 3 (unfortified industrial effluent) -- Statistical summary	52
Table C.4 -- Sample 4 (HPLC water) -- Statistical summary	53