

DIN EN ISO 4787:2022-03 (E)

Laboratory glass and plastic ware - Volumetric instruments - Methods for testing of capacity and for use (ISO 4787:2021)

Contents		Page
	European foreword	4
	Foreword	5
	Introduction	6
1	Scope	7
2	Normative references	7
3	Terms and definitions	8
4	Principle	8
5	Volume and reference temperature	8
	5.1 Unit of volume	8
	5.2 Reference temperature	8
6	Apparatus and calibration liquid	8
	6.1 Balance	8
	6.2 Measurement devices	8
	6.3 Calibration liquid	9
	6.4 Receiving vessel	9
7	Factors affecting the accuracy of volumetric instruments	9
	7.1 General	9
	7.2 Temperature	9
	7.2.1 Temperature of the volumetric instrument	9
	7.2.2 Temperature of calibration liquid	9
	7.3 Cleanliness of surface	10
	7.4 Conditions of used volumetric instruments	10
	7.5 Delivery time and waiting time	10
8	Setting the meniscus	11
	8.1 General	11
	8.2 Setting the meniscus	11
	8.2.1 Meniscus of transparent liquids	11
	8.2.2 Meniscus of opaque liquids	13
9	Calibration procedure	13
	9.1 General	13
	9.2 Test room	13
	9.3 Filling and delivery	13
	9.3.1 Volumetric flasks and measuring cylinders	13
	9.3.2 Pipettes adjusted to deliver	13
	9.3.3 Pipettes adjusted to contain	14
	9.3.4 Burettes adjusted to deliver	14
	9.3.5 Pycnometers	15
	9.4 Weighing	15
	9.5 Volume and uncertainty calculation	15
10	Procedure for use	16
	10.1 General	16
	10.2 Volumetric flasks (in accordance with ISO 1042 or ISO 5215)	17

10.3	Measuring cylinders (in accordance with ISO 4788 or ISO 6706)	17
10.4	Burettes (in accordance with ISO 385)	17
10.5	Pipettes	18
10.5.1	Pipettes adjusted to deliver (see ISO 648 and ISO 835, or other pipettes, e.g. plastic ones)	18
10.5.2	Pipettes adjusted to contain	18
10.6	Pycnometers	18
Annex A	(informative) Cleaning of volumetric glassware	19
Annex B	(informative) Cleaning of volumetric plasticware	20
Annex C	(normative) Calculation formulae and tables	21
Annex D	(informative) Coefficient of cubic thermal expansion	25
Annex E	(informative) Uncertainty estimation and repeatability calculation	26
Bibliography		27