

DIN ISO 16000-33:2017-12 (E)

Indoor air - Part 33: Determination of phthalates with gas chromatography/mass spectrometry (GC/MS) (ISO 16000-33:2017)

Contents		Page
National foreword		4
National Annex NA (informative) Bibliography		7
Foreword		10
Introduction		11
1 Scope		12
2 Normative references		12
3 Terms and definitions		12
4 Sampling methods and analytical apparatus		13
4.1 General		13
4.2 Sampling by adsorption with subsequent thermal desorption		13
4.2.1 Apparatus, operating materials and chemicals		13
4.2.2 Preparation of the thermal desorption tube		14
4.2.3 Sampling		14
4.3 Sampling by adsorption and subsequent solvent extraction		15
4.3.1 Apparatus, operating materials and chemicals		15
4.3.2 Preparation of Florisil® and the adsorption tubes		16
4.3.3 Hints to the application of Florisil®		17
4.3.4 Sampling		17
4.3.5 Sample conditioning		18
5 Calibration		18
5.1 General		18
5.2 Calibration of the thermal desorption method		19
5.3 Calibration of the solvent extraction method		19
6 Identification and quantification		19
6.1 Mass spectrometric analysis		19
7 Establishment of calibration curves and calculation of the analyte mass		25
7.1 Establishment of calibration curves		25
7.2 Calculation of the analyte mass		25
8 Calculation of indoor air concentrations		26
9 Performance characteristics		27
9.1 Detection limit		27
9.2 Quantification limit and problems related to the blank values		27
9.3 Reproducibility standard deviation and repeatability standard deviation		28
10 Quality assurance		30
10.1 Method verification and determination of blanks		30
10.1.1 Field blank value of the indoor air		30
10.1.2 Analytical laboratory blank value		30
10.2 Measures for blank value minimization		30
10.3 Documents		31
11 Interferences		31

Annex A (informative) General information on phthalates	32
Annex B (informative) Screening phthalates in solvent wipe tests	35
Annex C (informative) Screening phthalates in house dust	38
Annex D (informative) Practical example for the calibration of the thermal desorption method ...	42
Annex E (informative) Practical example for the calibration of the solvent extraction method using Florisil®	44
Annex F (informative) Practical example for the gas chromatography with thermal desorption ...	46
Annex G (informative) Practical example for the gas chromatography following solvent extraction	47
Annex H (informative) Problems related to the blank values	48
Annex I (informative) Example of a sampling protocol	49
Bibliography	50