

ISO 16773-2:2016-04 (E)

Electrochemical impedance spectroscopy (EIS) on coated and uncoated metallic specimens - Part 2: Collection of data

Contents		Page
Foreword		iv
1	Scope	1
2	Normative references	1
3	Terms and definitions	1
4	Principle	1
5	Electrochemical cell	2
5.1	General	2
5.2	Electrodes	3
5.3	Exposed area	3
5.4	Electrolyte	3
6	Procedure	4
6.1	Grounding	4
6.2	Shielding	4
6.3	Cell cable ground contacts	4
6.4	Local conditions	4
6.5	Measurement equipment characteristics	5
6.6	Confidence test	5
6.6.1	General	5
6.6.2	Use of interlaboratory test cells for confidence test	6
6.6.3	Error estimate and accuracy	6
6.7	Specimens	6
6.7.1	Preconditioning of specimens	6
6.7.2	Environmental control	6
6.8	Evaluation of laboratory and field coating specimens	6
6.9	Number of specimens and repeatability of results	7
7	Instrumental parameters	7
7.1	General	7
7.2	Conditioning potential and conditioning time	8
7.3	Frequency spectrum	8
7.4	Data point spacing and points per decade	8
7.5	D.C. potential control	8
7.6	Amplitude of perturbation	9
7.7	Current range settings	9
7.8	Data averaging settings	9
8	Data presentation	9
8.1	General	9
8.2	Bode plot	9
8.3	Nyquist plot	10
9	Exchange file format	11
10	Test report	11

Annex A (informative) Determination of maximum measurable impedance with the open- lead test	13
Annex B (normative)Dataexchangeformat	15
Annex C (informative) Considerations about the precise determination of the exposed area	19
Bibliography	24