

# ISO 25539-1:2017-02 (E)

## Cardiovascular implants - Endovascular devices - Part 1: Endovascular prostheses

---

<b>Contents</b>		<b>Page</b>
Foreword .....		v
Introduction .....		vi
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Normative references .....</b>	<b>2</b>
<b>3</b>	<b>Terms and definitions .....</b>	<b>2</b>
<b>4</b>	<b>General requirements for endovascular system .....</b>	<b>4</b>
4.1	Type of endovascular prosthesis .....	4
4.2	Materials and construction for endovascular system .....	4
4.3	Configuration and size designation for endovascular prosthesis .....	5
4.4	Intended clinical use for endovascular system .....	5
4.5	Balloon designation .....	6
<b>5</b>	<b>Intended performance .....</b>	<b>6</b>
<b>6</b>	<b>Design attributes .....</b>	<b>6</b>
6.1	General .....	6
6.2	Endovascular system .....	6
6.3	Endovascular prosthesis .....	6
6.4	Endovascular system and endovascular prosthesis .....	7
<b>7</b>	<b>Materials .....</b>	<b>7</b>
<b>8</b>	<b>Design evaluation .....</b>	<b>7</b>
8.1	General .....	7
8.2	Sampling .....	8
8.3	Conditioning of test samples .....	9
8.4	Reporting .....	9
8.5	Bench and analytical tests .....	10
8.5.1	Endovascular system and delivery system .....	10
8.5.2	Endovascular prosthesis .....	12
8.6	Preclinical in vivo evaluation .....	18
8.6.1	Purpose .....	18
8.6.2	Specific aims .....	18
8.6.3	Protocol considerations .....	19
8.6.4	Data acquisition .....	19
8.6.5	Test report and additional information .....	21
8.7	Clinical evaluation .....	21
8.7.1	Purpose .....	21
8.7.2	Specific aims .....	22
8.7.3	Protocol considerations .....	22
8.7.4	Data acquisition .....	23
8.7.5	Final report .....	26
<b>9</b>	<b>Post-market surveillance .....</b>	<b>27</b>
<b>10</b>	<b>Manufacturing .....</b>	<b>27</b>

<b>11</b>	<b>Sterilization .....</b>	<b>27</b>
<b>11.1</b>	<b>Products supplied sterile .....</b>	<b>27</b>
<b>11.2</b>	<b>Sterilization residuals .....</b>	<b>27</b>
<b>12</b>	<b>Packaging .....</b>	<b>28</b>
<b>12.1</b>	<b>Protection from damage in storage and transport .....</b>	<b>28</b>
<b>12.1.1</b>	<b>General .....</b>	<b>28</b>
<b>12.1.2</b>	<b>Unit container .....</b>	<b>28</b>
<b>12.1.3</b>	<b>Outer container .....</b>	<b>28</b>
<b>12.1.4</b>	<b>Shipping container .....</b>	<b>28</b>
<b>12.1.5</b>	<b>Maintenance of sterility in transit .....</b>	<b>28</b>
<b>12.2</b>	<b>Labelling .....</b>	<b>28</b>
<b>12.2.1</b>	<b>Container label .....</b>	<b>28</b>
<b>12.2.2</b>	<b>Record label .....</b>	<b>29</b>
<b>12.3</b>	<b>Instructions for use .....</b>	<b>29</b>
<b>12.3.1</b>	<b>General .....</b>	<b>29</b>
<b>12.3.2</b>	<b>Information and instructions for use for endovascular systems .....</b>	<b>29</b>
<b>Annex A (informative)</b>	<b>Relationship between testing requirements and device attributes and potential failure modes .....</b>	<b>31</b>
<b>Annex B (informative)</b>	<b>Description of clinical and device effects of failure .....</b>	<b>45</b>
<b>Annex C (informative)</b>	<b>Bench and analytical tests .....</b>	<b>49</b>
<b>Annex D (informative)</b>	<b>Test methods .....</b>	<b>57</b>
<b>Bibliography</b>	<b>.....</b>	<b>121</b>