

# DIN SPEC 58961:2013-04 (E)

## Technical Report - Thrombocyte diagnostics

---

<b>Contents</b>		<b>Page</b>
Foreword .....		3
1	Scope .....	4
2	Symbols and abbreviations .....	4
3	Pre-analytical procedures in platelet function testing .....	5
3.1	Patient preparation in platelet function testing .....	5
3.2	Blood extraction .....	5
3.2.1	Light transmission aggregometry (Born's method) (LTA) .....	5
3.2.2	Multiple electrode aggregometry (MEA) and impedance aggregometry .....	6
3.2.3	Platelet function analyser (PFA systems) .....	6
3.3	Determination of platelet count and packed cell volume .....	6
3.4	Transporting samples .....	6
3.5	Sample preparation .....	6
3.5.1	Light transmission aggregometry (Born's method) (LTA) .....	6
3.5.2	Impedance aggregometry .....	7
3.5.3	PFA systems .....	7
3.6	Reference range .....	7
4	Activation paths for platelet function testing .....	7
5	Platelet aggregation methods .....	9
5.1	Light transmission aggregometry (Born's method) (LTA) .....	9
5.1.1	General .....	9
5.1.2	Principle .....	9
5.1.3	Equipment, materials and reagents .....	9
5.1.4	Sample preparation .....	9
5.1.5	Procedure .....	9
5.1.6	Interpretation .....	10
5.1.7	Quality control .....	16
5.2	Aggregation methods -- Electrical impedance aggregometry (Cardinal's method) and multiple electrode aggregometry (MEA) .....	16
5.2.1	General .....	16
5.2.2	Principle .....	16
5.2.3	Equipment, materials and reagents .....	16
5.2.4	Sample preparation .....	17
5.2.5	Procedure .....	17
5.2.6	Interpretation .....	18
5.2.7	Quality control .....	19
6	The PFA system .....	20
6.1	General .....	20
6.2	Procedure .....	20
6.3	Interpretation of results .....	21
6.3.1	Expected results .....	21
6.3.2	Reference ranges and cut-off values .....	22
6.4	Clinical applications of PFA test cartridges .....	23
6.5	Quality control .....	23
7	Interfering factors .....	23
7.1	General .....	23

7.2	Method-specific interfering factors .....	24
7.2.1	LTA .....	24
7.2.2	Impedance aggregometry .....	24
7.2.3	PFA systems .....	24
	Bibliography .....	25