

DIN EN ISO 15883-7:2016-10 (E)

Washer-disinfectors - Part 7: Requirements and tests for washer-disinfectors employing chemical disinfection for non-invasive, non-critical thermolabile medical devices and healthcare equipment (ISO 15883-7:2016)

Contents

	Page
European foreword	4
Annex ZA (informative) Relationship between this European Standard and the essential requirements of Directive 93/42/EEC [OJ L 169] aimed to be covered	6
Foreword	13
Introduction	14
1 Scope	15
2 Normative references	15
3 Terms and definitions	16
4 Performance requirements	16
4.1 General	16
4.2 Cleaning	17
4.3 Disinfection	17
4.4 Final rinsing	19
4.5 Self-disinfection	19
4.6 Drying	20
4.7 Water treatment equipment	20
4.7.1 General	20
4.7.2 Disinfection of water treatment equipment	20
4.7.3 Maintenance of piping	21
5 Mechanical requirements	21
5.1 Materials — Design, manufacture, and assembly	21
5.2 Process verification	21
6 Testing for conformity	21
6.1 General	21
6.2 Test load	22
6.2.1 Loading with standard goods	22
6.2.2 Loading with special goods	22
6.3 Water used for rinsing following disinfection	22
6.4 Load dryness	22
6.4.1 General	22
6.4.2 Procedure	22
6.4.3 Results	22
6.5 Thermometric tests	22
6.5.1 General	22
6.5.2 Load temperature test	22
6.6 Chemical dosing tests	23
6.6.1 General	23
6.6.2 Reused process chemicals	23
6.7 Tests of cleaning efficacy	23
6.7.1 General	23
6.7.2 Materials	23
6.7.3 Procedure	24
6.7.4 Results	24

6.8	Test of disinfection efficacy	24
6.8.1	General.....	24
6.8.2	Preliminary tests on chemical disinfectants.....	24
6.8.3	Self-disinfection tests.....	25
6.8.4	Chemical disinfection of the load.....	26
7	Documentation.....	26
8	Information to be provided by the manufacturer.....	26
9	Marking, labelling, and packaging.....	27
10	Information to be requested from the purchaser by the manufacturer.....	27
Annex A	(normative) Summary of test programmes.....	28
Annex B	(normative) Methods for microbiological evaluation of disinfection of liquid transport system.....	30
Annex C	(normative) Tests for microbiological contamination of post-disinfection rinse water	35
Annex D	(normative) Preparation and evaluation of indicators for microbiological testing of the efficacy of chemical disinfection of the load.....	37
Annex E	(informative) Examples of test locations for the tests with biological indicators	40
Bibliography	44