

Contents

	Foreword	
	Introduction	
1	Scope	
2	Normative references	
3	Terms and definitions	
4	General considerations	
5	Outside the laboratory	
5.1	Specimen collection	
5.1.1	General	
5.1.2	Information about the specimen donor/patient	
5.1.3	Information about the specimen	
5.1.4	Specimen processing	
5.2	Fresh tissue transport requirements	
5.2.1	General	
5.2.2	Preparations for the transport	
5.2.3	During transport	
6	Inside the laboratory	
6.1	Information about the reception of the specimen	
6.2	Evaluation of the pathology of the specimen and selection of the sample(s)	
6.3	Freezing of the specimen or sample(s)	
6.4	Storage requirements	
6.5	Isolation of total protein	
6.5.1	General	
6.5.2	Using commercial kits	
6.5.3	Using the laboratories own protocols	
6.6	Quantity and quality assessment of isolated proteins	
6.7	Storage of isolated total protein	
Annex A	(informative) Quantitative protein examination demonstrates changes of protein amounts during cold ischemia ¹	Research by the EU FP7 SPIDIA project funded by the European Union Seventh Framework Programme [FP7/2007#2013] under grant agreement no 222916. For further information see www.spidia.eu .
A.1	Introduction	
A.2	Example	
A.2.1	General	
A.2.2	Experimental procedures	
A.2.2.1	General	
A.2.2.2	Tissues	
A.2.2.3	Protein analysis	
A.2.3	Results	
A.2.4	Further reading	