

Molecular in vitro diagnostic examinations — Specifications for pre-examination processes for venous whole blood — Part 2: Isolated genomic DNA

Contents

	Foreword	
	Introduction	
1	Scope	
2	Normative references	
3	Terms and definitions	
4	General considerations	
5	Outside the laboratory	
5.1	Specimen collection	
5.1.1	Information about the specimen donor/patient	
5.1.2	Selection of the venous whole blood collection tube by the laboratory	
5.1.3	Venous whole blood specimen collection from the donor/patient and stabilization procedures	
5.1.4	Information about the specimen and storage requirements at the blood collection facility	
5.1.4.1	General	
5.1.4.2	Using blood collection tubes with stabilizers	
5.1.4.3	Using blood collection tubes without stabilizers	
5.2	Transport requirements	
6	Inside the laboratory	
6.1	Specimen reception	
6.2	Storage requirements	
6.3	Isolation of the genomic DNA	
6.3.1	General	
6.3.2	Examination provider's instructions available	
6.3.3	Examination provider's instructions not available	
6.3.3.1	Using blood collection tubes without DNA stabilizers	
6.3.3.2	Using blood collection tubes with DNA stabilizers	
6.4	Quantity and quality assessment of isolated genomic DNA	
6.5	Storage of isolated genomic DNA	
6.5.1	General	
6.5.2	Genomic DNA isolated with commercially available kits	
6.5.3	Genomic DNA isolated with the laboratory's own protocols	
Annex A	(informative) Impact of pre-examination process steps on venous whole blood genomic DNA quality	
A.1	General information on operated experiments	
A.2	Influence of pre-examination variables (blood storage duration and temperature, and DNA isolation methods) on genomic DNA integrity	
A.3	Influence of genomic DNA integrity on an examination based on long PCR amplicons	
A.4	Influence of blood storage time on the genomic DNA integrity	
A.5	Influence of blood storage conditions on the performance of PCR tests based on short amplicons	