

# DIN 13222:2022-06 (E)

## Absorption capacity of absorbent incontinence aids until leakage - Test method for measuring the suction capacity by use of a test torso

---

<b>Contents</b>		<b>Page</b>
Foreword .....		4
Introduction .....		5
1 Scope .....		6
2 Normative references .....		6
3 Terms and definitions .....		7
4 Test method .....		9
4.1 Principle .....		9
4.2 Apparatus and assistive products .....		10
4.3 Conditioning .....		16
4.4 Procedure .....		16
4.4.1 Applying the test product on the mannequin .....		16
4.4.2 Pads (see Figure 11) .....		19
4.4.3 Belted pads (see Figure 14) .....		20
4.4.4 All-in-ones (see Figure 16) .....		21
4.4.5 Leakage detection .....		24
4.4.6 Performance of the procedure .....		24
4.4.7 Acceptance criteria .....		26
4.5 Test report .....		26
Annex A (informative) Validation of the test method .....		28
Annex B (informative) Example of a mounting frame .....		32
Annex C (informative) Procedure for calibrating the liquid stream .....		33
Bibliography .....		34
Figures Figure 1 -- Inclined position of the mannequin with male adapter .....		10
Figure 2 -- Mannequin dimensions (front and side view) for size 2 (medium, M) .....		11
Figure 3 -- Measuring position for determining the hardness .....		12
Figure 4 -- Dimension ranges of the male adapter .....		13
Figure 5 -- Assistive wedge .....		14
Figure 6 -- Micturition profile of the pump .....		16
Figure 7 -- Example of marking the assumed micturition point on various product versions .....		17
Figure 8 -- Centreline of the product aligned with the centreline between the legs of the mannequin .....		17
Figure 9 -- Product pouch in the crotch area .....		18

<b>Figure 10 -- Positioning the cuffs .....</b>	<b>19</b>
<b>Figure 11 -- Pad .....</b>	<b>19</b>
<b>Figure 12 -- The cuffs and elastic threads of the leg opening enclosing the male adapter .....</b>	<b>20</b>
<b>Figure 13 -- Product applied and held in place by mesh briefs .....</b>	<b>20</b>
<b>Figure 14 -- Belted pad .....</b>	<b>21</b>
<b>Figure 15 -- Applying belted pads .....</b>	<b>21</b>
<b>Figure 16 -- All-in-ones .....</b>	<b>21</b>
<b>Figure 17 -- Example fixation with closing strips for all-in-ones of product size M .....</b>	<b>22</b>
<b>Figure 18 -- Folding the side parts on a large all-in-one .....</b>	<b>23</b>
<b>Figure 19 -- Fixation of a large all-in-one on the mannequin .....</b>	<b>23</b>
<b>Figure 20 -- Leakage: test liquid exits the product .....</b>	<b>24</b>
<b>Figure 21 -- Positioning of the assistive wedge on the mounting frame (rear view) .....</b>	<b>25</b>
<b>Figure 22 -- Positioning of the assistive wedge in relation to the central axis of the mannequin .....</b>	<b>26</b>
<b>Figure B.1 -- Example of a mounting frame .....</b>	<b>32</b>
<b>Tables Table 1 -- Mannequin dimensions (front and side view) for size 2 (medium, M) .....</b>	<b>11</b>
<b>Table 2 -- Dimensions of the male adapter .....</b>	<b>12</b>
<b>Table 3 -- Flow data for the pump .....</b>	<b>15</b>
<b>Table A.1 -- Normal hip circumferences for all-in-ones within the specified size classes .....</b>	<b>28</b>
<b>Table A.2 -- Evaluation of the test method via an interlaboratory comparison of all-in-ones of size M and pads .....</b>	<b>28</b>
<b>Table A.3 -- Evaluation of the test method via an interlaboratory comparison of all-in-ones of sizes XS, S, L and XL .....</b>	<b>30</b>