

DIN 13222:2022-06 (E)

Absorption capacity of absorbent incontinence aids until leakage - Test method for measuring the suction capacity by use of a test torso

Contents		Page
Foreword		4
Introduction		5
1 Scope		6
2 Normative references		6
3 Terms and definitions		7
4 Test method		9
4.1 Principle		9
4.2 Apparatus and assistive products		10
4.3 Conditioning		16
4.4 Procedure		16
4.4.1 Applying the test product on the mannequin		16
4.4.2 Pads (see Figure 11)		19
4.4.3 Belted pads (see Figure 14)		20
4.4.4 All-in-ones (see Figure 16)		21
4.4.5 Leakage detection		24
4.4.6 Performance of the procedure		24
4.4.7 Acceptance criteria		26
4.5 Test report		26
Annex A (informative) Validation of the test method		28
Annex B (informative) Example of a mounting frame		32
Annex C (informative) Procedure for calibrating the liquid stream		33
Bibliography		34
Figures Figure 1 -- Inclined position of the mannequin with male adapter		10
Figure 2 -- Mannequin dimensions (front and side view) for size 2 (medium, M)		11
Figure 3 -- Measuring position for determining the hardness		12
Figure 4 -- Dimension ranges of the male adapter		13
Figure 5 -- Assistive wedge		14
Figure 6 -- Micturition profile of the pump		16
Figure 7 -- Example of marking the assumed micturition point on various product versions		17
Figure 8 -- Centreline of the product aligned with the centreline between the legs of the mannequin		17
Figure 9 -- Product pouch in the crotch area		18

Figure 10 -- Positioning the cuffs	19
Figure 11 -- Pad	19
Figure 12 -- The cuffs and elastic threads of the leg opening enclosing the male adapter	20
Figure 13 -- Product applied and held in place by mesh briefs	20
Figure 14 -- Belted pad	21
Figure 15 -- Applying belted pads	21
Figure 16 -- All-in-ones	21
Figure 17 -- Example fixation with closing strips for all-in-ones of product size M	22
Figure 18 -- Folding the side parts on a large all-in-one	23
Figure 19 -- Fixation of a large all-in-one on the mannequin	23
Figure 20 -- Leakage: test liquid exits the product	24
Figure 21 -- Positioning of the assistive wedge on the mounting frame (rear view)	25
Figure 22 -- Positioning of the assistive wedge in relation to the central axis of the mannequin	26
Figure B.1 -- Example of a mounting frame	32
Tables Table 1 -- Mannequin dimensions (front and side view) for size 2 (medium, M)	11
Table 2 -- Dimensions of the male adapter	12
Table 3 -- Flow data for the pump	15
Table A.1 -- Normal hip circumferences for all-in-ones within the specified size classes	28
Table A.2 -- Evaluation of the test method via an interlaboratory comparison of all-in-ones of size M and pads	28
Table A.3 -- Evaluation of the test method via an interlaboratory comparison of all-in-ones of sizes XS, S, L and XL	30