

# ISO 10282:2023-02 (E)

## Single-use sterile rubber surgical gloves - Specification

---

<b>Contents</b>		<b>Page</b>
Foreword .....		iv
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Normative references .....</b>	<b>1</b>
<b>3</b>	<b>Terms and definitions .....</b>	<b>1</b>
<b>4</b>	<b>Classification .....</b>	<b>2</b>
<b>4.1</b>	<b>General .....</b>	<b>2</b>
<b>4.2</b>	<b>Material .....</b>	<b>3</b>
<b>4.3</b>	<b>Design .....</b>	<b>3</b>
<b>4.4</b>	<b>Finish .....</b>	<b>3</b>
<b>5</b>	<b>Materials .....</b>	<b>3</b>
<b>6</b>	<b>Sampling and selection of test pieces .....</b>	<b>4</b>
<b>6.1</b>	<b>Sampling .....</b>	<b>4</b>
<b>6.2</b>	<b>Selection of test pieces .....</b>	<b>4</b>
<b>7</b>	<b>Requirements .....</b>	<b>4</b>
<b>7.1</b>	<b>Dimensions .....</b>	<b>4</b>
<b>7.2</b>	<b>Watertightness .....</b>	<b>5</b>
<b>7.3</b>	<b>Tensile properties .....</b>	<b>5</b>
<b>7.3.1</b>	<b>General .....</b>	<b>5</b>
<b>7.3.2</b>	<b>Force at break and elongation at break before accelerated ageing .....</b>	<b>5</b>
<b>7.3.3</b>	<b>Force at break and elongation at break after accelerated ageing .....</b>	<b>5</b>
<b>7.3.4</b>	<b>Force required to produce 300 % elongation .....</b>	<b>6</b>
<b>7.4</b>	<b>Sterility .....</b>	<b>6</b>
<b>8</b>	<b>Packaging .....</b>	<b>6</b>
<b>9</b>	<b>Marking .....</b>	<b>6</b>
<b>9.1</b>	<b>General .....</b>	<b>6</b>
<b>9.2</b>	<b>Inner package .....</b>	<b>6</b>
<b>9.3</b>	<b>Unit package .....</b>	<b>6</b>
<b>9.4</b>	<b>Multi-unit package .....</b>	<b>7</b>
<b>Annex A (normative) Measurement of dimensions (width, length and thickness) .....</b>		<b>8</b>
<b>Annex B (normative) Test method for watertightness .....</b>		<b>9</b>
<b>Bibliography .....</b>		<b>11</b>