

DIN EN 14180:2025-08 (E)

Sterilizers for medical purposes - Low temperature steam and formaldehyde sterilizers - Requirements and testing

| Contents | | Page |
|-------------------------|--|-------------|
| European foreword | | 6 |
| Introduction | | 8 |
| 1 | Scope | 9 |
| 2 | Normative references | 9 |
| 3 | Terms and definitions | 11 |
| 4 | General | 22 |
| 4.1 | Sterilizer type definition | 22 |
| 4.2 | Equipment development | 23 |
| 4.3 | Calibration | 24 |
| 5 | Equipment design and construction | 24 |
| 5.1 | Safety and security | 24 |
| 5.1.1 | General | 24 |
| 5.1.2 | Protective measures | 25 |
| 5.1.3 | Risk control and usability | 26 |
| 5.2 | Sterilizer chamber | 27 |
| 5.2.1 | Materials | 27 |
| 5.2.2 | Dimensions | 27 |
| 5.2.3 | Doors and interlocks of the sterilizer chamber | 28 |
| 5.2.4 | Temperature control and insulating the sterilizer chamber | 29 |
| 5.2.5 | Test connections | 29 |
| 5.3 | Further functional components | 31 |
| 5.3.1 | Pipework and fittings | 31 |
| 5.3.2 | Vaporizer | 31 |
| 5.3.3 | Evacuation system | 32 |
| 5.3.4 | Aeration system | 32 |
| 5.3.5 | Framework and panelling | 32 |
| 5.4 | Vibration | 33 |
| 5.5 | Ancillary equipment and components | 33 |
| 5.6 | Transport | 33 |
| 5.7 | User interfaces | 33 |
| 6 | Indicating, monitoring, controlling, and recording devices | 34 |
| 6.1 | Automatic controller | 34 |
| 6.2 | Automatic control | 35 |
| 6.3 | Control and monitoring system | 37 |
| 6.4 | Failures | 38 |
| 6.4.1 | General | 38 |
| 6.4.2 | Fault | 38 |
| 6.4.3 | Power failure | 39 |
| 6.4.4 | Other failures | 40 |
| 6.5 | Instrumentation | 40 |
| 6.5.1 | General | 40 |
| 6.5.2 | Temperature measuring devices | 41 |
| 6.5.3 | Pressure measuring devices | 42 |
| 6.5.4 | Time measuring devices | 43 |

| | | |
|--------|---|----|
| 6.5.5 | Sterilizing agent control and measuring devices | 43 |
| 6.6 | Indicating devices | 44 |
| 6.6.1 | General | 44 |
| 6.6.2 | Cycle parameter indicating devices | 44 |
| 6.6.3 | Cycle parameter indications | 46 |
| 6.6.4 | Status indicators and indications | 47 |
| 6.6.5 | Operating cycle counter | 48 |
| 6.7 | Recording systems | 48 |
| 6.7.1 | General | 48 |
| 6.7.2 | Records | 49 |
| 6.7.3 | Analogue presentation of records | 50 |
| 6.7.4 | Digital records | 50 |
| 6.7.5 | Data format for applied sterilant amount | 51 |
| 6.8 | Operating cycle | 51 |
| 6.8.1 | General | 51 |
| 6.8.2 | Leak test | 51 |
| 6.8.3 | Sterilization cycles | 52 |
| 7 | Services and local environment | 54 |
| 7.1 | General | 54 |
| 7.2 | Sterilizing agent and sterilant | 55 |
| 7.3 | Electrical supply | 56 |
| 7.4 | Water | 56 |
| 7.4.1 | Water used for sterilizer operation | 56 |
| 7.4.2 | Feed water | 56 |
| 7.5 | Steam | 56 |
| 7.6 | Vacuum | 56 |
| 7.7 | Drainage and discharges | 56 |
| 7.8 | Lighting | 57 |
| 7.9 | Compressed air | 57 |
| 7.10 | Air and inert gases | 57 |
| 7.11 | Ventilation and environment | 57 |
| 8 | Emissions | 58 |
| 8.1 | Electromagnetic emissions | 58 |
| 8.2 | Noise | 58 |
| 8.3 | Exhaust emissions | 59 |
| 8.4 | Heat emissions | 59 |
| 9 | Test instrumentation | 59 |
| 10 | Performance assessment | 59 |
| 10.1 | General | 59 |
| 10.2 | Attainment of conditions | 60 |
| 10.2.1 | Physical conditions | 60 |
| 10.2.2 | Sterilant and sterilizing agent | 62 |
| 10.2.3 | Microbicidal performance | 62 |
| 10.3 | Desorption efficacy | 62 |
| 10.4 | Load dryness | 62 |
| 10.5 | Pressure change | 62 |
| 11 | Information to be supplied | 63 |
| 11.1 | General | 63 |
| 11.2 | Information to be available prior to purchase | 63 |
| 11.3 | Information prior to installation | 64 |
| 11.4 | Marking and labelling | 65 |
| 11.5 | Instructions for use | 66 |
| 11.6 | Further information to be provided | 67 |
| 12 | Packaging | 68 |
| | Annex A (normative) Test methods | 69 |

| | | |
|---|--|------------|
| A.1 | General | 69 |
| A.2 | Test loads | 69 |
| A.3 | Test procedures | 71 |
| Annex B (normative) Sterilizer testing | | 78 |
| B.1 | General | 78 |
| B.2 | Type test | 78 |
| B.3 | Works test | 79 |
| B.4 | Installation qualification test provisions | 79 |
| B.5 | Operational qualification test provisions | 80 |
| Annex C (normative) Test equipment | | 81 |
| C.1 | Pressure instrumentation | 81 |
| C.2 | Temperature sensors | 81 |
| C.3 | Temperature recording instrument | 81 |
| C.4 | Pressure recording instrument | 82 |
| C.5 | Penetration type test device (PTTD) | 83 |
| C.6 | Residues challenge device | 83 |
| C.7 | Biological indicators and systems | 84 |
| C.8 | Sterile barrier systems | 84 |
| Annex D (normative) Determination of formaldehyde residuals in a filter indicator | | 85 |
| D.1 | Procedure for sample preparation | 85 |
| D.2 | Analysis of formaldehyde contents in filter indicator | 85 |
| Annex E (informative) Application of Formaldehyde in LTSF-processes | | 88 |
| E.1 | Formaldehyde concentrations, physico-chemical conditions during a cycle | 88 |
| E.2 | Formaldehyde residuals on medical devices | 90 |
| Annex F (informative) Illustrations of the interrelationship between control and recording | | 93 |
| F.1 | Introduction | 93 |
| F.2 | Illustration 1 | 95 |
| F.3 | Illustration 2 | 96 |
| F.4 | Illustration 3 | 97 |
| Annex G (informative) Additional information on protective measures | | 100 |
| Annex H (informative) Environmental aspects | | 101 |
| H.1 | Environmental aspects regarding the life cycle of LTSF sterilizers | 101 |
| H.2 | Formaldehyde (brief description) | 101 |
| H.3 | Environmental impact | 101 |
| Annex ZA (informative) Relationship between this European Standard and the General Safety and Performance Requirements of Regulation (EU) 2017/745 aimed to be covered | | 103 |
| Bibliography | | 125 |