

# DIN EN 16206:2012-05 (E)

**Animal feeding stuffs - Determination of arsenic by hydride generation atomic absorption spectrometry (HGAAS) after microwave pressure digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide)**

---

| <b>Contents</b>   |  | <b>Page</b> |
|---|--|-------------|
| Foreword .....  |  | 3           |
| 1   | Scope .....                            | 4           |
| 2   | Normative references .....             | 4           |
| 3   | Principle .....                        | 4           |
| 4   | Reagents .....                         | 5           |
| 5   | Apparatus and equipment .....          | 6           |
| 6   | Procedure .....                        | 7           |
| 6.1   | General .....                          | 7           |
| 6.2   | Preparation of the test solution ..... | 7           |
| 6.3   | Measurement of the test solution ..... | 8           |
| 7   | Calculation .....                      | 9           |
| 8   | Precision .....                        | 10          |
| 8.1   | Introduction .....                     | 10          |
| 8.2   | General .....                          | 10          |
| 8.3   | Repeatability .....                    | 10          |
| 8.4   | Reproducibility .....                  | 10          |
| 9   | Test report .....                      | 11          |
| Annex A (informative) Results of the inter-laboratory test .....  |  | 12          |
| Annex B (informative) Flowchart - Determination of arsenic by hydride generation atomic absorption spectrometry (HGAAS) after microwave pressure digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide) .....   |  | 13          |
| Annex C (informative) Alternative digestion procedure with the same digestion efficiency to ensure complete mineralization of all organic and inorganic arsenic species for HGAAS measurement: dry ashing with magnesium oxide and magnesium nitrate as ashing reagents ..... |  | 15          |
| Bibliography .....  |  | 16          |