

# DIN EN 16159:2012-04 (E)

**Animal feeding stuffs - Determination of selenium by hydride generation atomic absorption spectrometry (HGAAS) after microwave digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide)**

---

<b>Contents</b>	<b>Page</b>
Foreword .....	3
1 Scope .....	4
2 Normative references .....	4
3 Principle .....	4
4 Reagents .....	5
5 Apparatus and equipment .....	6
6 Procedure .....	7
6.1 General .....	7
6.2 Preparation of the test solution .....	7
6.3 Measurement of the test solution .....	8
6.3.1 Pre-dilution of the test solution .....	8
6.3.2 Pre-reduction of the (pre-diluted) test solution .....	8
6.3.3 Spectrometer settings of the flow-injection hydride generation atomic absorption spectrometer (HGAAS) .....	9
6.3.4 HGAAS determination .....	9
7 Calculation .....	10
8 Precision .....	10
8.1 Introduction .....	10
8.2 General .....	11
8.3 Repeatability .....	11
8.4 Reproducibility .....	11
9 Test report .....	11
Annex A (informative) Results of the inter-laboratory test .....	12
Annex B (informative) Flowchart - Determination of selenium by hydride generation atomic absorption spectrometry (HGAAS) after microwave digestion (digestion with 65 % nitric acid and 30 % hydrogen peroxide) .....	13
Annex C (informative) Alternative digestion procedure with the same digestion efficiency: Acid digestion with a mixture of 65 % nitric acid and 70 % perchloric acid (7:3 by volume) at atmospheric pressure .....	14
C.1 Preparation of selenium calibration solutions .....	14
C.2 Preparation of the test solution .....	14
C.3 Measurement of the test solution .....	15
C.4 Ten rules for automated wet ashing with perchloric acid .....	15
Bibliography .....	16