

# DIN EN 16274:2021-11 (E)

## Method for analysis of allergens - Quantification of an extended list of 57 suspected allergens in ready to inject fragrance materials by gas chromatography mass spectrometry

---

### Contents

	Page
European foreword .....	4
Introduction .....	5
1 Scope .....	6
2 Normative references .....	6
3 Terms and definitions .....	6
4 Principle .....	6
5 Reagents .....	7
5.1 Solvents .....	7
5.1.1 Methyl pivalate .....	7
5.1.2 Tert-Butyl Methyl Ether (MTBE) .....	7
5.2 Reference samples (suspected allergens) .....	7
5.2.1 General .....	7
6 Apparatus .....	11
6.1 Gas chromatograph equipped with flame ionization detector (GC-FID) .....	11
6.2 Gas chromatograph coupled to a mass spectrometer (GC-MS) .....	11
6.2.1 General .....	11
6.2.2 GC-MS system .....	11
6.3 Capillary columns for GC .....	11
6.4 Analytical balance .....	12
7 Stock and sample solutions - preparation and storage .....	12
7.1 General information .....	12
7.1.1 Choice of the solvent .....	12
7.1.2 Miscellaneous .....	13
7.2 Preparation of stock solutions from reference samples .....	13
7.3 Preparation of internal standard solutions and calibration solutions .....	13
7.3.1 General .....	13
7.3.2 Stock Solutions of Internal Standards (1g/kg equivalent) .....	13
7.3.3 Calibration Solutions Concentrations .....	13
7.3.4 Preparation of Calibration solutions .....	14
8 Mass spectrometer acquisition conditions .....	16
8.1 Establishment of retention times and SIM chromatograph windows .....	16
8.2 SIM Window - Set up and criteria .....	17
8.3 GC-MS scan mode verification .....	17
9 Sample analysis .....	17
9.1 General .....	17
9.2 Sample preparation for analysis .....	17
9.3 Sequence .....	18
9.3.1 Blank .....	18
9.3.2 Calibration solutions .....	18

<b>9.3.3</b>	<b>Check solution .....</b>	<b>18</b>
<b>9.3.4</b>	<b>Samples .....</b>	<b>19</b>
<b>10</b>	<b>Data validation and treatment .....</b>	<b>19</b>
10.1	General .....	19
10.2	Examination of Q values .....	19
10.3	Relative intensity of characteristic ions in SIM .....	20
10.4	Relative ion intensity using scan data .....	20
10.5	Data verification scheme and reporting of final concentration .....	21
10.6	GC-MS scan mode verification .....	21
10.7	Mass spectrum evaluation in SCAN mode .....	21
<b>11</b>	<b>Test report .....</b>	<b>21</b>
<b>Annex A (informative) Determination of calibrant and reference sample purity .....</b>		<b>22</b>
<b>Annex B (informative) GC capillary column parameters .....</b>		<b>29</b>
<b>Annex C (informative) SIM Ions for SIM or SIM-SCAN Use .....</b>		<b>30</b>
<b>Annex D (informative) Example of SIM windows .....</b>		<b>37</b>
<b>Annex E (informative) Example of chromatograms .....</b>		<b>46</b>
<b>Annex F (normative) Decisional tree for quantification of suspected allergens .....</b>		<b>48</b>
<b>Annex G (informative) Preparation of stock solutions from reference samples .....</b>		<b>49</b>
<b>Annex H (informative) Calibration methods and approach .....</b>		<b>54</b>
<b>Annex I (informative) Quantification of allergens .....</b>		<b>56</b>
<b>Annex J (informative) Suspected allergens selected as analytes targets - Rationale .....</b>		<b>59</b>
<b>Annex K (informative) Other general information .....</b>		<b>69</b>
<b>Bibliography .....</b>		<b>70</b>