

# DIN EN 16923:2023-01 (E)

## Foodstuffs - Determination of T-2 toxin and HT-2 toxin in cereals and cereal products for infants and young children by SPE clean up and HPLC-MS/MS

---

<b>Contents</b>		<b>Page</b>
European foreword .....		3
Introduction .....		4
1 Scope .....		5
2 Normative references .....		5
3 Terms and definitions .....		5
4 Principle .....		5
5 Reagents .....		5
6 Apparatus and equipment .....		8
7 Procedure .....		9
7.1 Preparation of the test sample .....		9
7.2 Preparation of the solid phase column .....		9
7.3 Extraction of T-2 toxin and HT-2 toxin .....		9
7.4 Clean-up by solid phase filtration .....		10
7.5 LC-MS/MS-analysis .....		10
7.6 Identification .....		10
8 Calculation .....		11
9 Precision .....		12
9.1 General .....		12
9.2 Repeatability .....		12
9.3 Reproducibility .....		12
10 Test report .....		13
Annex A (informative) Example chromatograms (API 4000TM) .....		14
Annex B (informative) Example conditions for suitable LC-MS/MS systems .....		18
B.1 System settings for SCIEX API 4000TM and SCIEX API 4000TM QTrap .....		18
B.1.1 Settings for chromatography .....		18
B.1.2 Detector parameters .....		18
B.2 System settings for SCIEX API 2000 .....		20
B.2.1 Settings for chromatography .....		20
B.2.2 Detector parameters .....		21
B.3 System settings for SCIEX API 3000TM .....		22
B.3.1 Settings for chromatography .....		22
B.3.2 Detector parameters .....		22
B.4 System settings for Micromass Quattro LC .....		23
B.4.1 Settings for chromatography .....		23
B.4.2 Detector parameters .....		24
Annex C (informative) Precision data .....		25
Bibliography .....		27