

# DIN EN ISO 15213-2:2024-05 (E)

## Microbiology of the food chain - Horizontal method for the detection and enumeration of *Clostridium* spp. - Part 2: Enumeration of *Clostridium perfringens* by colony-count technique (ISO 15213-2:2023)

---

<b>Contents</b>	<b>Page</b>
European foreword.....	3
Foreword.....	4
Introduction.....	6
<b>1 Scope</b> .....	<b>7</b>
<b>2 Normative references</b> .....	<b>8</b>
<b>3 Terms and definitions</b> .....	<b>8</b>
<b>4 Principle</b> .....	<b>9</b>
4.1 General.....	9
4.2 Preparation of dilutions.....	9
4.3 Enumeration.....	9
4.4 Confirmation.....	9
<b>5 Culture media and reagents</b> .....	<b>10</b>
<b>6 Equipment and consumables</b> .....	<b>10</b>
<b>7 Sampling</b> .....	<b>10</b>
<b>8 Preparation of test sample</b> .....	<b>11</b>
<b>9 Procedure</b> .....	<b>11</b>
9.1 General.....	11
9.2 Test portion, initial suspension and dilutions.....	11
9.3 Heat treatment to select spores.....	11
9.4 Inoculation and incubation.....	12
9.5 Enumeration of typical colonies.....	12
9.6 Confirmation of <i>C. perfringens</i> .....	13
9.6.1 Selection of colonies for confirmation.....	13
9.6.2 Acid phosphatase test.....	13
9.6.3 Sulfite indole motility (SIM) agar test.....	13
9.6.4 Differentiation between human pathogenic and non-pathogenic <i>C. perfringens</i> strains (optional).....	14
9.6.5 Interpretation.....	14
<b>10 Expression of results</b> .....	<b>14</b>
<b>11 Validation of the method</b> .....	<b>14</b>
11.1 Interlaboratory study.....	14
11.2 Performance characteristics.....	14
<b>12 Test report</b> .....	<b>15</b>
<b>13 Quality assurance</b> .....	<b>16</b>
<b>Annex A (normative) Flow diagram of the procedure</b> .....	<b>17</b>
<b>Annex B (normative) Culture media and reagents</b> .....	<b>19</b>
<b>Annex C (informative) Method validation studies and performance characteristics</b> .....	<b>24</b>
<b>Annex D (informative) Molecular differentiation between pathogenic and non-pathogenic     <i>C. perfringens</i></b> .....	<b>27</b>
<b>Bibliography</b> .....	<b>49</b>