

# DIN EN ISO 11816-2:2025-08 (E)

## Milk and milk products - Determination of alkaline phosphatase activity - Part 2: Fluorimetric method for cheese (ISO 11816-2:2024)

---

### Contents

Page

Forewords	iv
<b>1</b> Scope	<b>1</b>
<b>2</b> Normative references	<b>1</b>
<b>3</b> Terms and definitions	<b>1</b>
<b>4</b> Principle	<b>1</b>
<b>5</b> Reagents	<b>2</b>
<b>6</b> Apparatus	<b>3</b>
<b>7</b> Sampling	<b>4</b>
<b>8</b> Preparation of test sample	<b>4</b>
<b>9</b> Procedure	<b>5</b>
9.1 Verification of instrument performance	5
9.1.1 General	5
9.1.2 Daily instrument tests	5
9.1.3 Using FLM200	5
9.1.4 Using FLM300	5
9.1.5 Controls	6
9.2 Reagent controls to test the suitability of ready-to-use working substrate ( <a href="#">5.3</a> )	6
9.3 Calibration	6
9.3.1 General	6
9.3.2 Using FLM200	7
9.3.3 Using FLM300	7
9.4 Determination	7
9.5 Test-sample-related controls	8
9.5.1 Recommended negative and positive control tests	8
9.5.2 Interfering substance test	9
9.5.3 Heat-stable microbial ALP control test	9
<b>10</b> Calculation and expression of results	<b>9</b>
10.1 Calibration ratio	9
10.2 Calculation	10
10.2.1 Supernatant	10
10.2.2 Cheese	10
10.3 Expression of test results	10
<b>11</b> Precision	<b>11</b>
11.1 Interlaboratory study	11
11.2 Repeatability	11
11.3 Reproducibility	11
<b>12</b> Test report	<b>11</b>
<b>Annex A</b> (informative) <b>Interlaboratory study</b>	<b>12</b>
<b>Annex B</b> (informative) <b>Examples of preparation of a test sample</b>	<b>14</b>
<b>Bibliography</b>	<b>16</b>