

# ISO 15154:2003-04 (E)

## Tobacco - Determination of the content of reducing carbohydrates - Continuous-flow analysis method

---

<b>Contents</b>		<b>Page</b>
<b>1</b>	<b>Scope .....</b>	<b>1</b>
<b>2</b>	<b>Principle .....</b>	<b>1</b>
<b>3</b>	<b>Reagents .....</b>	<b>1</b>
<b>4</b>	<b>Apparatus .....</b>	<b>2</b>
<b>5</b>	<b>Procedure .....</b>	<b>3</b>
<b>6</b>	<b>Calculation .....</b>	<b>3</b>
<b>7</b>	<b>Repeatability and reproducibility .....</b>	<b>4</b>
<b>8</b>	<b>Test report .....</b>	<b>4</b>
<b>Annex A (informative) Prevention of precipitation of PAHBAH .....</b>		<b>5</b>
<b>Annex B (informative) Example of continuous-flow analyser .....</b>		<b>6</b>
<b>Bibliography .....</b>		<b>7</b>