

DIN EN 13946:2014-07 (E)

Water quality - Guidance for the routine sampling and preparation of benthic diatoms from rivers and lakes

| Contents | | Page |
|------------------------------|--|-------------|
| Foreword | | 3 |
| Introduction | | 4 |
| 1 | Scope | 5 |
| 2 | Terms and definitions | 5 |
| 3 | Principle | 6 |
| 4 | Equipment | 6 |
| 4.1 | Equipment for field sampling | 6 |
| 4.2 | Laboratory equipment | 6 |
| 5 | Reagents | 6 |
| 6 | Procedure | 8 |
| 6.1 | Choice of substratum | 8 |
| 6.2 | Sample site selection | 8 |
| 6.3 | Sampling methods | 8 |
| 6.3.1 | Moveable natural hard surfaces | 8 |
| 6.3.2 | Method for sampling vertical man-made surfaces in situ | 9 |
| 6.3.3 | Use of artificial substrata | 10 |
| 6.3.4 | Sample collection from submerged macrophytes and macroalgae | 10 |
| 6.3.5 | Sample collection from emergent macrophytes | 10 |
| 6.4 | Preparation prior to microscopic examination | 11 |
| 6.4.1 | Preservation and preliminary laboratory treatment | 11 |
| 6.4.2 | Methods for cleaning diatoms | 11 |
| 6.4.3 | Preparation of permanent slides | 11 |
| Annex A (informative) | Methods for cleaning diatoms for microscopic examination | 13 |
| A.1 | General | 13 |
| A.2 | Method 1: Hot hydrogen peroxide | 13 |
| A.2.1 | Apparatus | 13 |
| A.2.2 | Reagents | 13 |
| A.2.3 | Procedure | 13 |
| A.3 | Method 2: Cold hydrogen peroxide | 14 |
| A.3.1 | Apparatus and reagents | 14 |
| A.3.2 | Procedure | 14 |
| A.4 | Method 3: Hot hydrogen peroxide with potassium dichromate | 14 |
| A.4.1 | Apparatus | 14 |
| A.4.2 | Reagents | 14 |
| A.4.3 | Procedure | 14 |
| A.5 | Method 4: Cold acid (permanganate) method of cleaning | 15 |
| A.5.1 | Apparatus | 15 |
| A.5.2 | Reagents | 15 |
| A.5.3 | Procedure | 15 |
| Bibliography | | 17 |