

DIN EN ISO 19040-2:2023-12 (E)

Water quality - Determination of the estrogenic potential of water and waste water - Part 2: Yeast estrogen screen (A -YES, *Arxula adenivorans*) (ISO 19040-2:2018)

Contents	Page
European foreword	4
Foreword	5
1 Scope	6
2 Normative references	7
3 Terms and definitions	7
4 Principle	9
5 Interferences	9
6 Apparatus and materials	9
7 Reagents, media and test strains	10
8 Sampling and samples	14
8.1 General	14
8.2 Bottles and material for sampling	14
8.3 Bottles and material pre-cleaning	15
8.4 Sampling procedure	15
8.5 Transport of samples	15
8.6 Pretreatment of samples	15
8.7 Storage of samples	16
9 Procedure	16
9.1 Test set up	16
9.1.1 Preparation of the reference dilution series	16
9.1.2 Reactivation of the yeast	16
9.1.3 Negative control	17
9.1.4 Blank replicate	17
9.1.5 Sample dilution	17
9.1.6 Field blank	18
9.1.7 Plate set up	18
9.1.8 Inoculation of the test plate	18
9.2 Measurement	19
9.2.1 Measurement of the reporter gene activity	19
9.2.2 Measurement of the cell density	20
9.3 Calculations	20
9.3.1 Background correction	20
9.3.2 Calculation of the relative growth	21
9.3.3 Calculations for assessment of sample blanks	22
9.3.4 Calculation of the reporter gene induction	24
10 Validity criteria	27
11 Assessment criteria	28
12 Test report	28
13 Verification	28
Annex A (informative) Plate set up	30
Annex B (informative) Lyophilization of <i>Arxula adenivorans</i> cells	31
Annex C (informative) Scheme of test principle	33
Annex D (informative) Test set up for chemicals and extracts	34
Annex E (informative) Preparation of dilution series	35

Annex F (informative) Performance data	36
Annex G (informative) Statistical assessment	46
Annex H (informative) Calculation of estradiol equivalents	47
Annex I (informative) Alternative test design for EEQ determination	50
Annex J (informative) Measurement of the lowest ineffective dilution (LID) of a waste water — A simplified evaluation for testing of waste water	51
Annex K (informative) Example for statistical evaluation	53
Bibliography	60