

DIN EN 14105:2011-07 (E)

Fat and oil derivatives - Fatty Acid Methyl Esters (FAME) - Determination of free and total glycerol and mono-, di-, triglyceride contents

Contents	Page
Foreword	3
1 Scope	4
2 Principle	4
3 Reagents	4
4 Apparatus	5
5 Preparation of solutions	6
5.1 1,2,4-Butanetriol stock solution, 1 mg/ml	6
5.2 Glycerol stock solution, 0,5 mg/ml	6
5.3 Standard glycerides stock solution, 2,5 mg/ml	6
5.4 Commercial mixture of monoglycerides	6
5.5 Calibration solutions	6
6 Sampling	7
7 Procedure	7
7.1 Operating conditions	7
7.2 Analysis of the calibration solutions	7
7.3 Analysis of the commercial mixture of monoglycerides	7
7.4 Preparation and analysis of the samples	8
7.5 Identification	8
7.6 Calibration	8
7.7 Column performance control	8
8 Determination of results	9
8.1 Integration of the peaks	9
8.2 Glycerol calibration function	9
8.3 Free glycerol	9
8.4 Glycerides	10
8.5 Total glycerol	10
9 Expression of results	10
10 Precision	10
10.1 Interlaboratory test	10
10.2 Repeatability	11
10.3 Reproducibility	11
11 Test report	11
Annex A (informative) Sample chromatogram	12
Annex B (normative) Calibration function calculation	16
Annex C (informative) Worked example	18
Annex D (informative) Results of the interlaboratory trial	19
Bibliography	21