

DIN EN 14105:2011-07 (E)

Fat and oil derivatives - Fatty Acid Methyl Esters (FAME) - Determination of free and total glycerol and mono-, di-, triglyceride contents

Contents		Page
Foreword		3
1	Scope	4
2	Principle	4
3	Reagents	4
4	Apparatus	5
5	Preparation of solutions	6
5.1	1,2,4-Butanetriol stock solution, 1 mg/ml	6
5.2	Glycerol stock solution, 0,5 mg/ml	6
5.3	Standard glycerides stock solution, 2,5 mg/ml	6
5.4	Commercial mixture of monoglycerides	6
5.5	Calibration solutions	6
6	Sampling	7
7	Procedure	7
7.1	Operating conditions	7
7.2	Analysis of the calibration solutions	7
7.3	Analysis of the commercial mixture of monoglycerides	7
7.4	Preparation and analysis of the samples	8
7.5	Identification	8
7.6	Calibration	8
7.7	Column performance control	8
8	Determination of results	9
8.1	Integration of the peaks	9
8.2	Glycerol calibration function	9
8.3	Free glycerol	9
8.4	Glycerides	10
8.5	Total glycerol	10
9	Expression of results	10
10	Precision	10
10.1	Interlaboratory test	10
10.2	Repeatability	11
10.3	Reproducibility	11
11	Test report	11
Annex A (informative) Sample chromatogram		12
Annex B (normative) Calibration function calculation		16
Annex C (informative) Worked example		18
Annex D (informative) Results of the interlaboratory trial		19
Bibliography		21