

DIN EN 17126:2019-02 (E)

Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of sporicidal activity of chemical disinfectants in the medical area - Test method and requirements (phase 2, step 1)

Contents

Page

European foreword	4
Introduction	5
1 Scope	6
2 Normative references	6
3 Terms and definitions	6
4 Requirements	6
5 Test method	7
5.1 Principle	7
5.2 Materials and reagents	8
5.2.1 Test organisms	8
5.2.2 Culture media and reagents	8
5.3 Apparatus and glassware	14
5.3.1 General	14
5.3.2 Usual microbiological laboratory equipment	14
5.4 Preparation of test organism suspensions and product test solutions	16
5.4.1 Test organism suspensions (test and validation suspension)	16
5.4.2 Product test solutions	21
5.5 Procedure for assessing the sporicidal activity of the product	21
5.5.1 General	21
5.5.2 Dilution-neutralization method	23
5.5.3 Membrane filtration method	25
5.5.4 Modified method for ready-to-use products	27
5.6 Experimental data and calculation	29
5.6.1 Explanation of terms and abbreviations	29
5.6.2 Calculation	30
5.7 Verification of methodology	34
5.7.1 General	34
5.7.2 Control of weighted mean counts	34
5.7.3 Basic limits	34
5.8 Expression of results and precision	35
5.8.1 Reduction	35
5.8.2 Control of active and non-active product test solution (5.4.2)	35
5.8.3 Limiting test organism and sporicidal concentration	35
5.8.4 Precision, repetitions	36
5.9 Interpretation of results - conclusion	36
5.9.1 General	36
5.9.2 Sporicidal activity for surface disinfection products	36
5.9.3 Sporicidal activity for instrument disinfection products	36
5.9.4 Sporicidal activity for textile disinfection products	36
5.9.5 Qualification for certain fields of application	36
5.10 Test report	36
Annex A (informative) Referenced strains in national collections	39

Annex B (informative) Neutralizers and Rinsing Liquids	40
Annex C (informative) Graphical representation of test procedures.....	42
C.1 Dilution-neutralization method	42
C.2 Membrane filtration method	44
C.3 Dilution-neutralization method (modified method for ready-to-use products)	46
C.4 Membrane filtration method (modified method for ready-to-use products)	48
Annex D (informative) Example of a typical test report	50
Annex E (informative) Precision of the test result	54
Annex F (informative) Graphical representation of spore preparation	55
Annex G (informative) Guide for preparation of reference test solution	56
Annex H (informative) Example for a titration Peracetic acids and hydrogen peroxide	57
H.1 General	57
H.2 Equipment	57
H.3 Procedure	57
Bibliography	59